Stat Profile Prime Plus® VET Calibrator Cartridge With Creatinine

Stat Profile Prime Plus® VET Καλιβριεστική με Κρεατίνη, Φυσίγγιο μέσου βαθμονόμησης VET με κρεατινίνη Stat Profile Prime Plus®, Cartucho de calibrador con creatinina Stat Profile Prime Plus® VET, Cartouche de calibratageur avec créatinine Stat Profile Prime Plus® VET, Cartucho de calibrador Stat Profile Prime Plus® VET com creatinina, Stat Profile Prime Plus® VET kalibrator patron kreatininnel, Stat Profile Prime Plus® VET קריאטינין כויל עם קריאטינין, Stat Profile Prime Plus® VET キャリブレータ (クレアチニン), 크레아티닌 사용 Stat Profile Prime Plus® VET 보정기 카트리지, Stat Profile Prime Plus® 含肌酐的校准试剂盒

LOT 24257071

2026-03-04

CAL A	≥ 495 mL	CAL B	≥ 250 mL	CAL C	≥ 245 mL	CAL D	≥ 250 mL	CAL F	≥ 850 mL	SOLN R	≥ 500 mL
HCO ₃	11 mmol/L	pH	6.834 @ 37°C	HCO ₃	22 mmol/L	pH	7.354 @ 37°C	PO ₂	103.8 mmHg	KCl	2 mol/L
Na ⁺	135.5 mmol/L	H ⁺	146.55 nmol/L	Cl ⁻	69.0 mmol/L	H ⁺	44.26 nmol/L	PO ₂	13.81 kPa	WASTE CONTAINER W	
K ⁺	4.0 mmol/L	Na ⁺	77.0 mmol/L	Glu	200 mg/dL						
Cl ⁻	108.0 mmol/L	K ⁺	10.0 mmol/L	Glu	11.10 mmol/L						
iCa	1.04 mmol/L	iCa	1.92 mmol/L	Lac	10.0 mmol/L						
iCa	4.16 mg/dL	iCa	7.68 mg/dL	Lac	89.0 mg/dL						
iMg	0.50 mmol/L	iMg	1.50 mmol/L	BUN	50 mg/dL						
iMg	1.22 mg/dL	iMg	3.63 mg/dL	BUN	17.9 mmol/L						
Glu	80 mg/dL	tHb	15.4 g/dL	Creat	5.0 mg/dL						
Glu	4.44 mmol/L	tHb	154 g/L	Creat	0.44 mmol/L						
Lac	2.0 mmol/L	tHb	9.6 mmol/L	Creat	442 μmol/L						
Lac	17.8 mg/dL										
BUN	10 mg/dL										
BUN	3.6 mmol/L										
Creat	1.0 mg/dL										
Creat	0.09 mmol/L										
Creat	88.4 μmol/L										

EN
Description
 The Cartridge contains the following aqueous solutions within individual bags: 5 Calibrators (Calibrator A, B, C, D, and F), Reference Solution, and a flexible waste bag. Each bag includes a filament with a septum. The exposed bag filaments are arranged in a line along the rear of the cartridge box. The septa are pierced during the insertion of a line along the analyzer. A charge kit consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation.

Intended Use
 For the quantitative determination of Na⁺, K⁺, Cl⁻, iCa, iMg, BUN (Urea), Glucose, Lactate, Creatinine, pH, PCO₂, PO₂, SO₂, Hct, Hb, O₂Hb, COHb, MetHb, HtHb, and tHb in animal whole blood using the Stat Profile Prime Plus VET Analyzer.

Methodology
 Refer to Stat Profile Prime Plus VET Analyzer Instructions For Use Manual for Methodologies and Principles.

Composition
 Calibrator and solution volumes are listed in the table. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-T2).

Warnings and Cautions
 Refer to the Stat Profile Prime Plus VET Analyzer Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory reagents. Discard Cartridge in accordance with local guidelines. Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

Storage
 Store at 15 - 30°C. Expiration date is printed on the Calibrator Cartridge.

Directions For Use
Activation of Cartridge
 Cal A and Cal C must be charged with creatinine syringe prior to installation of the Calibrator Cartridge. The two syringes are labeled and color coded to correspond to the filaments on the carton for Cal A and Cal C.
 1. Hold the syringe with lip side down and remove protective cap.
 2. Attach one of the enclosed needle assemblies to the syringe. Remove the protective cover from the needle.
 3. Match the color and label of the syringe to the appropriate filament and insert needle.
 4. Slowly depress syringe plunger until the contents are dispensed.
 DO NOT PULL BACK ON THE PLUNGER TO FLUSH CONTENTS OF SYRINGE.
 5. Remove needle/syringe assembly from filament and discard in an appropriate sharps container.
 6. Repeat Steps 1-5 for the next Cal.
 7. Mix Cartridge well by gently shaking for 1 minute.

Refer to the Stat Profile Prime Plus VET Analyzer Instructions for Use Manual for complete directions. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use.

Traceability of Standards
 Analytes traced to NIST Standard Reference Materials wherever possible.

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UK

Description
The cartridge contains the following aqueous solutions within individual bags: 5 calibrators (calibrators C, D and F), reference solution and a flexible waste bag. Each bag includes an activator with a septum. The activator bags are exposed in a line on the back of the cartridge box. The septa are pierced during the insertion of the cartridge into the charging device consisting of 2 syringes filled with creatinine solution, for the activation of the cartridge before installation.

For indication
For the quantitative determination of Na⁺, K⁺, Cl⁻, iCa, iMg, BUN (urea), glucose, lactate, creatinine, pH, PCO₂, PO₂, SO₂, % Hct, iHb, O₂Hb, COHb, MetHb, HbH and iBil in whole animal blood using the Stat Profile Prime Plus VET analyser.

Methodology
For methodology and principles of analysis, see the manual of instructions for use of the Stat Profile Prime Plus VET analyser.

Composition
The volumes of calibrators and solutions are indicated in the table. It does not contain any human-derived substances; however, good laboratory practice should be followed when handling these materials. (REF. NCCLS DOCUMENT M29-12).

Warnings and precautions
For detailed instructions, see the manual of instructions for use of the Stat Profile Prime Plus VET analyser. Follow standard laboratory practices for handling reagents. Dispose of the cartridge in accordance with local regulations.

Storage
Store at 15-30°C. The expiry date is printed on the calibrator cartridge.

Instructions for use
Cartridge activation
Before installing the calibrator cartridge, load the calibrators A and C with the creatinine syringe. Dispose of the calibrator cartridge in accordance with local regulations.

1. Hold the syringe with the needle tip pointing downwards and remove the protective cap.
2. Connect one of the included water sets to the syringe. Remove the protective cap of the water.
3. Match the colour and label of the syringe with the appropriate fitting and insert the needle.
4. Push the plunger of the syringe slowly until the contents have been dispensed.
5. Remove the needle/syringe assembly from the CONTENT OF THE SYRINGE.
6. Repeat steps 1-5 for the next Cal.
7. Mix the cartridge well by gently shaking it for 1 minute.

For detailed instructions, see the manual of instructions for use of the Stat Profile Prime Plus VET analyser. The cartridge is ready for use.

Conformance with standards
Results are traceable to NIST reference materials where possible.

US

Description
The cartridge contains the following aqueous solutions within individual bags: 5 calibrators (calibrators C, D and F), reference solution and a flexible waste bag. Each bag includes an activator with a septum. The activator bags are exposed in a line on the back of the cartridge box. The septa are pierced during the insertion of the cartridge into the charging device consisting of 2 syringes filled with creatinine solution, for the activation of the cartridge before installation.

For indication
For the quantitative determination of Na⁺, K⁺, Cl⁻, iCa, iMg, BUN (urea), glucose, lactate, creatinine, pH, PCO₂, PO₂, SO₂, % Hct, iHb, O₂Hb, COHb, MetHb, HbH and iBil in whole animal blood using the Stat Profile Prime Plus VET analyser.

Methodology
For methodology and principles of analysis, see the manual of instructions for use of the Stat Profile Prime Plus VET analyser.

Composition
The volumes of calibrators and solutions are indicated in the table. It does not contain any human-derived substances; however, good laboratory practice should be followed when handling these materials. (REF. NCCLS DOCUMENT M29-12).

Warnings and precautions
For detailed instructions, see the manual of instructions for use of the Stat Profile Prime Plus VET analyser. Follow standard laboratory practices for handling reagents. Dispose of the cartridge in accordance with local regulations.

Storage
Store at 15-30°C. The expiry date is printed on the calibrator cartridge.

Instructions for use
Cartridge activation
Before installing the calibrator cartridge, load the calibrators A and C with the creatinine syringe. Dispose of the calibrator cartridge in accordance with local regulations.

1. Hold the syringe with the needle tip pointing downwards and remove the protective cap.
2. Connect one of the included water sets to the syringe. Remove the protective cap of the water.
3. Match the colour and label of the syringe with the appropriate fitting and insert the needle.
4. Push the plunger of the syringe slowly until the contents have been dispensed.
5. Remove the needle/syringe assembly from the CONTENT OF THE SYRINGE.
6. Repeat steps 1-5 for the next Cal.
7. Mix the cartridge well by gently shaking it for 1 minute.

For detailed instructions, see the manual of instructions for use of the Stat Profile Prime Plus VET analyser. The cartridge is ready for use.

Conformance with standards
Results are traceable to NIST reference materials where possible.

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