REF 57831, 57832, 57833, 57834, 57835, 65142

Rx Only IVD

15°C -30°C CE

# Stat Profile Prime Plus® Calibrator Cartridge With Creatinine

Stat Profile Prime Plus® Kalibrierkassette mit Kreatinin, Stat Profile Prime Plus® Φυσιγγίο μέσου βαθμονόμησης με κρεατινίνη, Cartucho de calibrador con creatinina Stat Profile Prime Plus®, Cartouche de calibrateurs Stat Profile Prime Plus® avec créatinine, Cartuccia calibratore con creatinina Stat Profile Prime Plus®, Cartucho de calibrador Stat Profile Prime Plus® com creatinina, Stat Profile Prime Plus® kalibrátor patron kreatininnel, ιν्यां क्रैटिनीन तरीकों से मापने वाला Stat Profile Prime Plus®, Stat Profile Prime Plus® キャリブレーターカートリッジ (クリアーチン), 크레아티닌 사용 Stat Profile Prime Plus® 보정기 카트리지, Stat Profile Prime Plus® 含肌酐的校准试剂盒

LOT

24257071

2026-03-04

CAL A	CAL B	CAL C	CAL D	CAL F	SOLN R
HCO <sub>3</sub> 11 mmol/L	pH 6.834 @ 37°C	HCO <sub>3</sub> 22 mmol/L	pH 7.354 @ 37°C	PO <sub>2</sub> 103.8 mmHg	KCl 2 mol/L
Na <sup>+</sup> 135.5 mmol/L	H <sup>+</sup> 146.55 nmol/L	Cl <sup>-</sup> 69.0 mmol/L	H <sup>+</sup> 44.26 nmol/L	PO <sub>2</sub> 13.81 kPa	
K <sup>+</sup> 4.0 mmol/L	Na <sup>+</sup> 77.0 mmol/L	Glu 200 mg/dL			
Cl <sup>-</sup> 108.0 mmol/L	K <sup>+</sup> 10.0 mmol/L	Glu 11.10 mmol/L			
iCa 1.04 mmol/L	iCa 1.92 mmol/L	Lac 10.0 mmol/L			
iCa 4.16 mg/dL	iCa 7.68 mg/dL	Lac 89.0 mg/dL			
iMg 0.50 mmol/L	iMg 1.50 mmol/L	BUN 50 mg/dL			
iMg 1.22 mg/dL	iMg 3.63 mg/dL	BUN 17.9 mmol/L			
Glu 80 mg/dL	tHb 15.4 g/dL	Creat 5.0 mg/dL			
Glu 4.44 mmol/L	tHb 154 g/L	Creat 0.44 mmol/L			
Lac 2.0 mmol/L	tHb 9.6 mmol/L	Creat 442 µmol/L			
Lac 17.8 mg/dL					
BUN 10 mg/dL					
BUN 3.6 mmol/L					
Creat 1.0 mg/dL					
Creat 0.09 mmol/L					
Creat 88.4 µmol/L					

EN

**Description**  
The Cartridge contains the following aqueous solutions within individual bags: 5 Calibrators (Calibrator A, B, C, D, F), Reference Solution, and a flexible waste bag. Each bag includes a filament with a septum. The exposed bag filaments are arranged in a line along the rest of the cartridge box. The septa are pierced during the insertion of the cartridge into the analyzer. A charge kit consisting of 2 syringes filled with creatinine solution is included for activation prior to installation.

#### Intended Use

For the quantitative determination of Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, Ca<sup>2+</sup>, iMg, BUN (Urea), Glucose, Lactate, Creatinine, pH, PCO<sub>2</sub>, PO<sub>2</sub>, SO<sub>4</sub><sup>2-</sup>, Hct, tHb, HbP<sup>-</sup>, O<sub>2</sub>Hb, COHb, MetHb, HrHb, and tBil<sup>+</sup> in human blood using the Stat Profile Prime Plus Analyzer.

#### Methodology

Refer to Stat Profile Prime Plus Analyzer Instructions for Use Manual for Methodologies and Principles.

#### Compatibility

The calibrators are aqueous salt solutions containing analytes at concentrations listed in the table above. Calibrator A/B also contains dye. Certain non constituents of human origin, often good laboratory practice should be followed during the handling of these materials. (REF. NCCLS DOCUMENT M29-T2).

#### Warnings and Cautions

Intended for *in vitro* diagnostic use.  
Refer to the Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory reagents. Discard Cartridge in accordance with local guidelines.

Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

#### Storage

Store at 15 - 30 °C. Expiration date is printed on the Calibrator Cartridge.

#### Directions For Use

##### Activation of Cartridge

Cal A and Cal C must be charged with creatinine syringe prior to installation of the Calibrator Cartridges. The cartridges are labeled with color codes to correspond to the filaments on the cartridge for Cal A and Cal C.

- Hold the syringe with tip side down and remove protective cap.
- Attach one of the enclosed needle assemblies to the syringe. Remove the protective cover from the needle.
- Slowly depress the plunger until the contents are dispensed.
- DO NOT FULL BACK ON THE PLUNGER TO FLUSH CONTENTS OF SYRINGE.
- Remove needle/syringe assembly from filament and discard in an appropriate sharps container.
- Repeat Steps 1-5 for the next Cal.
- Mix cartridge well by gently shaking for 1 minute.

Refer to the Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use..

\*Not available in the USA or for Point-of-Care/Near-Patient testing use.

#### Traceability of Standards

Analyses traced to NIST Standard Reference Materials wherever possible.

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