Stat Profile Prime Plus® Blood Gas, CO-Oximeter, Chemistry Controls Auto-Cartridge with Creatinine

Cartucho automático con creatinina para controles de química, gases en sangre y cooxímetro Stat Profile Prime Plus®. Cartouche automatique de contrôles de gaz du sang/CO-oxymètre, chimie aéreins αιμάτος, CO-Oξύμετρο Stat Profile Prime Plus®, Cartuccia con creatinina per controlli automatici chimici per gas ematici/CO-ossimetría Stat Profile Prime Plus®, Cartucho automático de controles de química, de CO-oxímetro e de Gás no sangue Stat Profile Prime Plus® com creatinina, Stat Profile Prime Plus® vérág, CO-oximéter, kémiai kontrollok automatikus patron kreatininnel, de CO-oxímetro e de Gás no sangue Stat Profile Prime Plus® Co-Oximeter כרמיין, Stat Profile Prime Plus® 血液ガス、CO オキシメーター、生化学検査用コントロール自動カートリッジ(creatinnin), 케레이티난 사용 Stat Profile Prime Plus® 혈액 가스, CO-산소 농도계, 화학 조절제 자동 카트리지, Stat Profile Prime Plus® 血气、一氧化碳血氧仪、化学对照溶液自动试剂盒(含肌酐)

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EN

Product Description:
Analyze the cartridge material for monitoring the performance of pH, PO₂, PO₃, SO₂, hematocrit (Hct), hemoglobin (Hgb), total hemoglobin (HtB), total bilirubin (BilB), cytochrome c oxidase (CO₂), carboxyhemoglobin (COHb), metahemoglobin (MetHb) and deoxyhemoglobin (DyHb) in Levels 1, 2 and 3 as well as Control Level 4 and 5. A charge kit consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation. For use on UNICENTER™ with the analyzers Stat Profile Prime PLUS.

Intended Use:Designed for *in vitro* diagnostic use by healthcare professionals for monitoring the performance of the Stat Profile Prime Plus

Stat Profile Prime Analyzer. Units are FDA 510(k) cleared.

Methodology: Refer to Stat Profile Prime Plus Analyzer Instructions For Use Manual for Methodology and Principles.**Compositos:** Control Levels 1, 2 and 3 are buffered bicarbonate solutions containing dyes, salts and preservatives. Each level has a known pH and contains a known O₂, CO₂ and N₂ value. Control Level 4 and 5 are buffered solutions containing known concentrations of Na⁺, K⁺, Cl⁻, Ca²⁺, Mg²⁺, Glucose, Lactate, BUN (urea), Creatinine and preservatives. Each pouch contains a minimum of 100 mL. Components are of animal or human origin, however good laboratory practices should be followed during handling of these materials. (REF: NCCLS Document M29-T2).**Warnings and Cautions:**

DO NOT FREEZE. Mix the cartridge by gently inverting for several seconds. DO NOT SHAKE CARTRIDGE. Refer to Stat Profile Prime Analyzer Instructions for Use Manual for complete information.

Intended for *in vitro* diagnostic use. Follow standard practices for handling laboratory reagents.**Storage:** Store at 2-8°C (37-45°F). DO NOT FREEZE.**Directions for use:** Ensure the cartridge room temperature prior to installation.

Level 4 and Level 5 must be charged with creatinine syringes prior to installation of the Control Cartridge. The 2 syringes are labeled and coded to match the labels on the Control Cartridge.

Activation steps are as follows:

1. Hold the syringe with tip down and remove the protective cap.

2. Attach the needle to the end of the syringe and remove the protective cover from the needle.

3. Slowly depress syringe plunger until the contents are dispensed. DO NOT PUSH BACK ON THE PLUNGER TO FLUSH CARTRIDGE. Hold the syringe upright and insert the needle.

4. Remove needle/syringe assembly from the cartridge and discard in an appropriate sharps container.

5. Repeat Steps 1-5 for the next Control.

6. Mix the cartridge by gently inverting for 1 minute. Cartridge is ready for use.

Verify that the lot number on the Control Cartridge corresponds to the lot number on the cartridge. Refer to Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions.

Limitations: PO₂ values may vary inversely with temperature (approximately 1% /°C). Therefore, it is critical to follow the temperature guidelines described in "Directions for use". The expected Range values are specific for instruments and controls manufactured by Nova Biomedical. Once installed, each Stat Profile Prime Plus cartridge may be used for a maximum of 21 days from the initial installation date. If the cartridge is used after 21 days, the system will indicate the cartridge is invalid. Each cartridge may be inserted and removed from the analyzer a maximum of 3 times.**Traceability of Standards:** Total Hemoglobin (HtB) and Metahemoglobin (MetHb) are traceable by using Cyanmethemoglobin method. Carboxyhemoglobin (COHb) and Cytochrome c oxidase (CO₂) are traceable using Spectrophotometry. Analyses are traceable to NIST Standard Reference Materials.**Reference Intervals:** Control levels are formulated at normal and abnormal expected values in patient blood. The expected clinical range of these values in patient blood is referenced in Tietz, NW ed. 1986 Textbook of Clinical Chemistry, W.B. Saunders Co. Users may wish to determine Mean Values and Expected Ranges in their own laboratory.²**Expected Ranges:** The expected range for each parameter was determined at Nova Biomedical using replicate determinations on Nova analyzers. The expected range indicates the maximum deviation from the Mean Value that may be expected under differing laboratory conditions for instruments operating within specifications. Refer to Expected Ranges Table.¹Not available in the USA or for Point-of-Care/Patient Testing use.

NCCLS Document M29-T2.

2How to Define and Determine Reference Intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13.

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Waltham, MA 02454 U.S.A.