StatProfile Prime Ampule Control ES

CONTROL 1 2

LOT 24255049

2026-08-29

Expected Ranges, Erwartungsbereiche, Αναμενόμενα εύρη, Intervalos esperados,
Plages attendues, Intervalli previsti, Várható tartományok, 予測範囲,
예상 범위, Gamas previstas, 预期范围

LOT		CONTROL 1		CONTROL 2	
CONTROL 1	CONTROL 2	min - \bar{x} - max		min - \bar{x} - max	
24242036	24242037	2026-08-29		2026-08-29	
pH*		7.34 - 7.39 - 7.44		7.50 - 7.55 - 7.60	
H+*	nmol/L	45.7 - 40.7 - 36.3		31.6 - 28.2 - 25.1	
Hct	%	48 - 51 - 54		63 - 66 - 69	
Na+	mmol/L	140.2 - 145.2 - 150.2		114.1 - 117.6 - 121.1	
K+	mmol/L	3.45 - 3.65 - 3.85		5.98 - 6.28 - 6.58	
Cl-	mmol/L	101.1 - 105.1 - 109.1		80.6 - 84.6 - 88.6	
iCa	mg/dL	4.17 - 4.81 - 5.45		6.81 - 7.62 - 8.42	
iCa	mmol/L	1.04 - 1.20 - 1.36		1.70 - 1.90 - 2.10	
iMg	mg/dL	0.80 - 1.12 - 1.43		2.62 - 3.06 - 3.50	
iMg	mmol/L	0.33 - 0.46 - 0.59		1.08 - 1.26 - 1.44	

- * Only Available on the Stat Profile ES Comp Plus Analyzer
- * Disponible únicamente en el analizador Stat Profile ES Comp Plus
- * Nur am Stat Profile ES Comp Plus Analysegerät verfügbar
- * Διατίθεται μόνο στον Αναλυτή Stat Profile ES Comp Plus.
- * Disponible uniquement sur l'analyseur Stat Profile ES Comp Plus.
- * Disponibile esclusivamente sull'analizzatore Stat Profile ES Comp Plus
- * Csak a Stat Profile ES Comp Plus analizátoron áll rendelkezésre
- * Stat Profile ES Comp Plus 분석기에서만 사용 가능함
- * Apenas disponível no Stat Profile ES Comp Plus Analyzer
- * 仅供Stat Profile ES Comp Plus 分析器使用
- * Stat Profile ES Comp Plus分析装置でのみ使用できます

EC REP Nova Biomedical GmbH Hessenring 13A, Geb. G 64546 Mörfelden-Walldorf Germany

EN

Product Description
Contains assayed quality control material formulated at three levels for monitoring the measurement of pH, Hct, Na⁺, K⁺, Cl⁻, Li⁺, Ca²⁺, and Mg²⁺.

Intended Use
For *in vitro* diagnostic use to monitor the performance of the Stat Profile Prime ES, Stat Profile Prime ES Plus, and Stat Profile Prime ES Comp Plus Analyzers.**

Methodology
Refer to StatProfile Prime Analyzer Instruction for Use Manual for Methodology and Principles.

Composition
Nova Control is a pale straw colored liquid assayed quality control material. Contains organic and inorganic constituents in a bovine albumin base. A preservative has been added to the product; however it does not contain stabilizers, viscosity adjusters, drug analogs or other materials which might affect the performance of sensors. Controls contain a buffered solution with known concentrations of pH, Hct, Na⁺, K⁺, Cl⁻, Li⁺, Ca²⁺, and Mg²⁺. Each ampule contains a minimum volume of 1.2 mL. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials.¹

Warning and Precautions
For *in vitro* diagnostic use. Store at 2-8°C. Do not Freeze. Refer to the Nova Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory reagents.

Storage
Store at 2-8°C. Do not freeze. Each ampule has a lot number and expiration date printed on the label.

Directions for Use
Re-suspend the ampule contents. Nova Controls contain proteins, which may settle with storage. Gently rotate each ampule between fingers or hands to mix contents prior to use. Avoid shaking, excess foaming may occur. Snap open the ampule (protecting fingers with gauze or gloves), aspirate solution with sample probe or transfer solution into sample cups. Verify that the lot number appearing on the expected ranges table is the same as the lot number printed on the material.

Limitations
Assigned assay values will be adversely affected by improper handling or storage, or if the integrity of the product packaging has been compromised prior to use. Cloudiness or particulate matter in the solution is an indication of possible deterioration. For use on Nova Biomedical Analyzers Only. Performance characteristics have not been established for use on other manufacturer's analyzers.

This product is warranted to perform as described in this product insert, and Nova Biomedical disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall Nova Biomedical be liable for any consequential damages arising out of the aforesaid expressed warranty.

Traceability of Standards
Analytes are traced to NIST Standard reference materials.

¹NCCLS Document M29-T2.
²The Stat Profile Prime ES Plus and Stat Profile Prime ES Comp Plus Analyzers are not currently CE marked.

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Waltham, MA 02454 U.S.A.