Stat Profile Prime Plus® Blood Gas, CO-Oximeter, Chemistry Controls Auto-Cartridge with Creatinine

Cartucho automático con creatinina para controles de química, gases en sangre y cooxímetro Stat Profile Prime Plus®. Cartouche automatique de contrôles de gaz du sang/CO-oxymètre, chimie Stat Profile Prime Plus® avec créatine. Stat Profile Prime Plus® Auto-Kassette mit Kreatinin für Blutgas-, CO-Oximeter- und Blutchemiekontrollen. Αυτόματο φουίγγιο υαίκοιο ελέγχου χημείας αερίων αίματος, CO-Οξόμετρο Stat Profile Prime Plus®, Cartuccia con creatinina per controlli automatici chimici per gas ematici/CO-ossimetria Stat Profile Prime Plus®, Cartucho automático de controles de química, de CO-oxímetro e de Gás no sangue Stat Profile Prime Plus® com creatinina, Stat Profile Prime Plus® vérgáz, CO-oxíméter, kémiai kontrollok automatikus patron kreatininnal, קריאטינין, סטאט פרופיל פרייμ פלוס סטאט פרופיל פרייμ פלוס קו-אוקסיμέטר, Stat Profile Prime Plus® 血液ガス, CO オキシメーター, 生化学検査用コントロール自動カートリッジ(クレアチニン), 크레아티닌 사용 Stat Profile Prime Plus® 혈액 가스, CO-산소 농도계, 화학 조절제 자동 카트리지, Stat Profile Prime Plus® 血气、一氧化二碳血氣儀、化学对照溶液自动试剂盒(含肌酐)

LOT 24292044

2026-04-01

CONTROL 1 2 3 4 5

Table with 6 columns: Parameter, CONTROL 1, CONTROL 2, CONTROL 3, CONTROL 4, CONTROL 5. Rows include pH, H+, PCO2, PO2, SO2, Hct, Na+, K+, Cl-, iCa, cCa, iMg, Glu, Hb, HbA1c, Urea, Creatinine, HbF, tHb, iHb, O2Hb, COHb, MetHb, HHb, tBil, iBil.

EN Product Description: NovaBio quality control material for monitoring the performance of pH, PCO2, PO2, SO2, hematocrit (Hct), fetal hemoglobin (HbF), total hemoglobin (Hb), total bilirubin (tBil), urea/hemoglobin (CO-Hb), methemoglobin (MetHb), and deoxyhemoglobin (HHb) in Levels 1, 2, and 3 as well as Na+, K+, Cl-, Ca, Mg, Glucose, Lactate, BUN (Urea), and Creatinine in Level 4 and 5. A change in constituents of 2 ranges from one creatinine solution is included for control activation prior to installation. For use with Stat Profile Prime Plus Analyzers ONLY.

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EN Directions for use: Ensure controls are room temperature prior to installation. Level 4 and Level 5 must be changed with creatinine syringes prior to installation of the Control Cartridge. The 2 syringes are labeled and color coded to correspond to the fillers on the cartridge. Activate the cartridge as follows: 1. Insert the syringe with the side down and remove protective cap. 2. Attach one of the enclosed needle assemblies to the syringe. Remove the protective cover from the needle. 3. Match the color and label of the syringe to the appropriate filler and insert the syringe. 4. Engage the plunger to the cartridge. DO NOT PUSH BACK ON THE PLUNGER TO FLUSH CONTENTS OF SYRINGE. 5. Remove needle/syringe assembly from filler and discard in an appropriate sharps container. 6. Repeat Steps 1-6 for the next Control. 7. Mix Cartridge well by gently inverting for 1 minute. Cartridge is ready for use. Verify that the Lot Number on the Expected Ranges Table corresponds to the Lot Number in the cartridge. Refer to Stat Profile Prime Plus Analyzers Instructions for Use Manual for complete directions.

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EN *Not available in the USA or for Point-of-Care/Point-of-Use Testing use. (NCCLS Document M29-12) Show to Define and Determine Reference Intervals in the clinical laboratory, approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13

ES *No disponible en EE. UU. o para uso en puntos de atención del paciente. (NCCLS Document M29-12) Show to Define and Determine Reference Intervals in the clinical laboratory, approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13

FR *Non disponible aux États-Unis ou pour les points d'intervention/diagnostic délocalisés. (NCCLS Document M29-12) Show to Define and Determine Reference Intervals in the clinical laboratory, approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13

NOVA biomedical logo and address: Nova Biomedical GmbH, Hessenring 13A, Geb. G 64546 Mörfelden-Walldorf Germany, 200 Prospect Street, Waltham, MA 02454 U.S.A.