NOVA Linearity Creatinine/BUN/Hct Levels 1-4

Solución de linealidad para creatinina/BUN/Hct de Nova, Niveles 1 a 4, Niveaux de linéarité 1-4 NOVA pour créatinine/BUN/Hct, NOVA Linearitätslösung für Kreatinin/BUN/Hct, Level 1-4, Επίπεδα 1-4 Κρεατινίνης Γραμμικότητας/BUN/Hct NOVA, Soluzioni di linearità NOVA per creatinina/BUN/Hct Livelli 1-4, Linearidade da NOVA Níveis 1-4 de creatinina/BUN/Hct, NOVA linearitás kreatinin/BUN/Hct 1-4. szint, 1-4Hct/BUN NOVA Linearity Creatinine רמת, NOVA リニアリティークレアチニン/BUN/Hct レベル 1-4, NOVA 선형성 크레아티닌/BUN/Hct 레벨 1-4, NOVA线性肌酐/尿素氮/红细胞比容1-4

LOT 24298035

2025-10-22

		Expected Ranges, Rangos esperados, Plages attendues, Erwartungsbereiche, Αναμενόμενο Εύρος, Intervalli previsti, Intervalli previstos, Várt tartományok, Expected Ranges (HE), צפויים טווחים, 予測範囲, 예상 범위, 预期范围值			
		LEVEL 1	LEVEL 2	LEVEL 3	LEVEL 4
LOT		24296026	24296037	24297023	24297024
		2025-10-22	2025-10-22	2025-10-23	2025-10-23
BUN	mg/dL	59 - 71 - 83	38 - 48 - 58	17 - 24 - 31	2 - 5 - 8
BUN	mmol/L	21.1 - 25.3 - 29.6	13.6 - 17.1 - 20.7	6.1 - 8.6 - 11.1	0.7 - 1.8 - 2.9
Urea	mg/dL	126.6 - 152.3 - 178.0	81.5 - 103.0 - 124.4	36.5 - 51.5 - 66.5	4.3 - 10.7 - 17.2
Urea	mmol/L	21.1 - 25.3 - 29.6	13.6 - 17.1 - 20.7	6.1 - 8.6 - 11.1	0.7 - 1.8 - 2.9
Creatinine	mg/dL	7.2 - 9.7 - 12.2	5.5 - 7.0 - 8.5	3.4 - 4.6 - 5.8	0.5 - 1.0 - 1.5
Creatinine	mmol/L	0.64 - 0.86 - 1.08	0.49 - 0.62 - 0.75	0.30 - 0.41 - 0.51	0.04 - 0.09 - 0.13
Creatinine	μmol/L	640 - 860 - 1080	490 - 620 - 750	300 - 410 - 510	40 - 90 - 130
Hct	%	12 - 15 - 18	26 - 29 - 32	44 - 47 - 50	62 - 65 - 68

EN

Product Description

Aqueous solutions containing Na⁺, K⁺, Cl⁻, HCO₃⁻, Urea, and Creatinine intended for monitoring the performance of the Creatinine, BUN and hematocrit channels.

Intended Use

For *in vitro* diagnostic use with Stat Profile Prime Plus Analyzers to verify calibration, analytical linearity, estimate test imprecision and detect systematic analytical deviations that may arise from calibrator or analytical instrument variation.

Methodology

Refer to Nova Analyzer Instructions for Use for Methodology and Principles of the testing procedure.

Composition

Buffered solutions containing Na⁺, K⁺, Cl⁻, HCO₃⁻, Urea, and Creatinine. Each ampule contains a 1.7 mL volume. The concentrations of the active ingredients have been adjusted to give conductivity characteristics equivalent to a known hematocrit value in whole blood. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-T2).

Warnings and Cautions:

Intended for *in vitro* diagnostic use. Follow standard practices for handling laboratory reagents. Discard in accordance with local requirements.

Storage

Store at -15°C. After thawing, the unopened solution is stable for up to 14 days when stored refrigerated at 2-8°C.

Directions for use

Prepare the solution by thawing to room temperature. Before opening, shake the ampule for 10 seconds. Snap open the ampule (protecting fingers with gloves or gauze). Linearity Standards are recommended within this multipack. It is recommended that each laboratory establish the upper and lower limits of its reportable range for each sample type and devise a policy for processing specimens that fall outside of this range.

Limitations

For use on Nova Biomedical Analyzers only. Performance characteristics not established for use on other manufacturer's analyzers.

Traceability of Standards

Analytes are traced to NIST Standard Reference Materials.

Expected Ranges

The expected range indicates the maximum deviations from the mean value that may be expected under differing laboratory conditions for instruments operating within specifications. The expected ranges for the analytes listed in the table above were verified using replicate determinations on Nova analyzers.

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