Nova Primary Ampuled Controls 1 and 2

Contrôles en ampoule Nova Primary 1 et 2, Controlli 1 e 2 in fiala Nova Primary, Controles en ampolla 1 y 2 de Nova Primary, Controles 1 e 2 em ampolas Nova Primary, Nova Primary Kontrolllösungen 1 und 2 in Ampullen

CONTROL 1 2

LOT 23355017 2025-07-08

LOT		Expected Ranges, Plages attendues, Intervalli previsti, Rangos esperados, Intervalos previstos, Erwartungsbereiche	
CONTROL 1	24008043	2025-07-08	
CONTROL 2	24008042	2025-07-08	
		CONTROL 1	CONTROL 2
		min - \bar{x} - max	min - \bar{x} - max
Glu	mg/dL	68 - 72 - 76	180 - 192 - 204
Glu	mmol/L	3.8 - 4.0 - 4.2	10.0 - 10.7 - 11.3

EN

Product Description

Aqueous quality control material for monitoring the performance of Glucose for use with the Nova Primary Glucose Analyzer. It is formulated at 2 levels.

Intended Use

Intended for *in vitro* diagnostic use by healthcare professionals use for monitoring the performance of the Nova Primary Glucose Analyzer.

Methodology

Refer to Nova Primary Glucose Analyzer Instruction for Use Manual for Methodology and Principles.

Composition

Controls are buffered salt solutions containing known concentrations of Glucose and preservatives. Each ampule contains a minimum volume of 1.7 mL. Controls contain no constituents of human origin, however, good laboratory practices should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-T21').

Warning and Precautions

Intended for *in vitro* diagnostic use. DO NOT FREEZE. Follow standard practices for handling laboratory reagents. Discard the unused portion in accordance with local guidelines.

Storage

Store at 15-30°C. DO NOT FREEZE. Each ampule has a Lot Number and Expiration Date printed on the label.

Directions for Use

Shake ampule before opening, snap open ampule (protecting fingers with gauze or glove). Discard the unused portion in accordance with local guidelines. Verify that the Lot Number on the Expected Ranges Table corresponds to the Lot Number on the ampule. Refer to the Nova Primary Glucose Analyzer Instructions for Use Manual for complete instructions.

Limitations

The Expected Range values are specific for instruments manufactured by Nova Biomedical.

Traceability of Standards

Analytes are traced to NIST Standard reference materials.

Reference Intervals

Concentrations are formulated at normal and abnormal expected values in patient blood. The expected clinical range of these values in patient blood is referenced in Tietz, NW ed. 1986 Textbook of Clinical Chemistry, W.B. Saunders Co. Users may wish to determine Mean Values and Expected Ranges in their laboratory.²

Expected Ranges

The expected range was determined at Nova Biomedical using replicate determinations on Nova analyzers. The expected range indicates the maximum deviations from the Mean Value that may be expected under differing laboratory conditions for instruments operating within specifications. Refer to Expected Ranges Table.

1NCCLS Document M29-T2.

2How to Define and Determine Reference Intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13.

EC REP Nova Biomedical GmbH Hessenring 13A, Geb. G 64546 Mörfelden-Walldorf Germany

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