REF 57831, 57832, 57833, 57834, 57835, 65142

Rx Only IVD 15°C 30°C CE

# Stat Profile Prime Plus® Calibrator Cartridge With Creatinine

Stat Profile Prime Plus® Kalibrirkassette mit Kreatinin, Stat Profile Prime Plus® Φυσίγγιο μέσου βαθμονόμησης με κρεατινίνη, Cartucho de calibrator con creatinina Stat Profile Prime Plus®, Cartouche de calibrateurs Stat Profile Prime Plus® avec créatinine, Cartuccia calibratore con creatinina Stat Profile Prime Plus®, Cartucho de calibrator Stat Profile Prime Plus® com creatinina, Stat Profile Prime Plus® kalibrátor patron kreatininnel, מוסגית הכיול עם קריאטינין Stat Profile Prime Plus®, Stat Profile Prime Plus® キャリブレーターカートリッジ (クレアチニン), 크레아티닌 사용 Stat Profile Prime Plus® 보정기 카트리지, Stat Profile Prime Plus® 含肌酐的校准试剂盒, Stat Profile Prime Plus® Kreatinin İçeren Calibrator Cartridge

LOT 24331049

2026-03-18

CAL A		CAL B		CAL C		CAL D		CAL F		SOLN R	
HCO <sub>3</sub>	11 mmol/L	pH	6.834 @ 37°C	HCO <sub>3</sub>	22 mmol/L	pH	7.354 @ 37°C	PO <sub>2</sub>	102.9 mmHg	KCl	2 mol/L
Na <sup>+</sup>	135.5 mmol/L	H <sup>+</sup>	146.55 nmol/L	Cl <sup>-</sup>	69.0 mmol/L	H <sup>+</sup>	44.26 nmol/L	PO <sub>2</sub>	13.69 kPa		<b>WASTE CONTAINER</b> <b>W</b>
K <sup>+</sup>	4.0 mmol/L	Na <sup>+</sup>	77.0 mmol/L	Glu	200 mg/dL						
Cl <sup>-</sup>	108.0 mmol/L	K <sup>+</sup>	10.0 mmol/L	Glu	11.10 mmol/L						
iCa	1.04 mmol/L	iCa	1.92 mmol/L	Lac	10.0 mmol/L						
iCa	4.16 mg/dL	iCa	7.68 mg/dL	Lac	89.0 mg/dL						
iMg	0.50 mmol/L	iMg	1.50 mmol/L	BUN	50 mg/dL						
iMg	1.22 mg/dL	iMg	3.63 mg/dL	BUN	17.9 mmol/L						
Glu	80 mg/dL	tHb	15.4 g/dL	Creat	5.0 mg/dL						
Glu	4.44 mmol/L	tHb	154 g/L	Creat	0.44 mmol/L						
Lac	2.0 mmol/L	tHb	9.6 mmol/L	Creat	442 µmol/L						
Lac	17.8 mg/dL										
BUN	10 mg/dL										
BUN	3.6 mmol/L										
Creat	1.0 mg/dL										
Creat	0.09 mmol/L										
Creat	88.4 µmol/L										

### EN Description

The Cartridge contains the following aqueous solutions within individual bags: 5 Calibrators (Calibrator A, B, C, D, and F), Reference Solution, and a flexible waste bag. Each bag includes a fliant with a septum. The exposed bag fliants are arranged in a line along the rear of the cartridge box. The septa are pierced during the insertion of the cartridge into the analyzer. A charge kit consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation.

**Intended Use**  
For the quantitative determination of Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, iCa, iMg, BUN (Urea), Glucose, Lactate, Creatinine, pH, PCO<sub>2</sub>, PO<sub>2</sub>, SO<sub>2</sub>%, Hct, Hb, HbF\*, O<sub>2</sub>Hb, COHb, Methb, HbB, and tHb\* in human blood using the Stat Profile Prime Plus Analyzer.

**Methodology**  
Refer to Stat Profile Prime Plus Analyzer Instructions For Use Manual for Methodologies and Principles.

**Composition**  
The calibrators are aqueous salt solutions containing analytes at concentrations listed in the table above. Calibrator B also contains dye. Contain no constituents of human origin, however, good laboratory practice should be followed during the handling of these materials. (REF. NOCLS DOCUMENT M29-72).

**Warnings and Cautions**  
Intended for in vitro diagnostic use. Refer to the Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory reagents. Discard Cartridge in accordance with local guidelines.

**Storage**  
Store at 15 - 30°C. Expiration date is printed on the Calibrator Cartridge.

**Directions For Use**  
**Activation of Cartridge**  
Cal A and Cal C must be charged with creatinine syringe prior to installation of the Calibrator Cartridge. The two syringes are labeled and color coded to correspond to the fittings on the carton for Cal A and Cal C.

- Hold the syringe with tip side down and remove protective cap.
- Attach one of the enclosed needle assemblies to the syringe. Remove the protective cover from the needle.
- Match the color and label of the syringe to the appropriate fliant and insert needle.
- Slowly depress syringe plunger until the contents are dispensed.
- DO NOT PULL BACK ON THE PLUNGER TO FLUSH CONTENTS OF SYRINGE.
- Remove needle/syringe assembly from fliant and discard in an appropriate sharps container.
- Repeat Steps 1-5 for the next Cal.
- Mix Cartridge well by gently shaking for 1 minute.

Refer to the Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use.

\*Not available in the USA or for Point-of-Care/Near-Patient Testing use.

**Traceability of Standards**  
Analyses traced to NIST Standard Reference Materials wherever possible.

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