REF 48559, 41549

IVD 15°C - 30°C CE

Stat Profile® pHox Ultra/Critical Care Xpress CO-Oximeter Calibrator Cartridge with Bilirubin and Deproteinizing Solution

Kalibrator-Patrone mit Bilirubin und Deproteinisierungslösung, Πακέτο βαθμονομητή με χοληρυθρίνη και διάλυμα αποπρωτεϊνισμού, Cartucho de calibrador con bilirubina y solución desproteinizante, Cartouche de calibreur avec bilirubine et solution de déproteinisation, Cartuccia calibratore con bilirubina e soluzione desproteinizzante, Cartucho de calibrador com solução de desproteinização e bilirubina, ビリルビンおよび除蛋白液付きキャリブレーター・カートリッジ, Calibrator Cartridge de bilirubina 및 단백질 제거액, 含有胆红素和去蛋白溶液的定标品盒

LOT 24291037 2025-10-17

Table with 5 columns: tHb CALIBRATOR ≥70mL, DIL ≥180mL, SOLN F ≥620mL, Deproteinizing Solution ≥90mL, WASTE CONTAINER W. Each column contains detailed instructions in multiple languages.

Description For use in calibrating the CO-Oximeter on the Stat Profile pHox Ultra/Critical Care Xpress (COx) Analyzer. The above table lists the concentrations of the constituents for the calibrator, flush, deproteinizing and diluent solutions contained in the Calibrator Cartridge.

Intended Use For the quantitative determination of total hemoglobin, oxyhemoglobin, carboxyhemoglobin, methemoglobin, deoxyhemoglobin and total bilirubin in human blood using the Stat Profile pHox Ultra/Critical Care Xpress CO-Oximeter Analyzer.

Methodology Refer to the Stat Profile pHox Ultra/Critical Care Xpress Analyzer Instructions for Use Manual for Test Methodologies and Principles.

Composition Solution volumes are listed in the table above. Contains no constituents of human origin, however good laboratory practice should be followed during handling of these materials.

Warnings and Cautions Intended for in vitro diagnostic use. Refer to the Stat Profile pHox Ultra/Critical Care Xpress Analyzer Instructions for Use Manual for complete directions. Cartridge is ready for use.

Directions for Use Refer to the Stat Profile pHox Ultra/Critical Care Xpress Analyzer Instructions for Use Manual for complete directions for installation of the cartridge. Cartridge is ready for use.

Limitations For use on Nova Biomedical Analyzers only. If not stored at room temperature and out of direct sunlight, the solution may deteriorate.

Traceability of Standards Traceable to precise hemoglobin measurement using the reference cyanmethemoglobin method.

Reference ICSH, Recommendations for Hemoglobinometry in Human Blood. British Journal of Haematology, 13 (suppl) pp 71-75, 1967.

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EC REP Nova Biomedical GmbH Hessenring 13A, Geb. G 64546 Mörfelden-Walldorf Germany
200 Prospect Street Waltham, MA 02454 U.S.A. www.novabiomedical.com