Stat Profile Prime Plus® VET Blood Gas, CO-Oximeter, Chemistry Controls Auto-Cartridge with Creatinine

Cartucho automático con creatinina para controles de química, gases en sangre y cooxímetro Stat Profile Prime Plus® VET, Cartouche automatique de contrôles de gaz du sang/CO-oxymètre, chimie Stat Profile Prime Plus® VET avec créatine, Stat Profile Prime Plus® VET Auto-Kassette mit Kreatinin für Blutgas-, CO-Oximeter- und Blutchemiekontrollen, Αυτόματο φιατίγγιο με κρεατινίνη επιπέδου ελέγχου χημείας αερίων αίματος, CO-Οξυμέτρο Stat Profile Prime Plus® VET, Cartuccia con creatinina per controlli automatici chimici per gas ematici/CO-ossimetria Stat Profile Prime Plus® VET, Cartucho automático de controles de química, de CO-oxímetro e de gás no sangue Stat Profile Prime Plus® VET com creatinina, Stat Profile Prime Plus® VET vérsgáz, CO-oxíméter, kémiai kontrollok automatikus patron kreatininnel, (מסוימת בקרינה אולטראסונית עם קריאטינין) Stat Profile Prime Plus® VET Chemistry Controls Auto-Cartridge, Stat Profile Prime Plus® VET 血液ガス, CO オキシメーター, 生化学検査用コントロール自動カートリッジ(クレアチニン), 크레아티닌 사용 Stat Profile Prime Plus® VET 혈액 가스, CO-산소 농도계, 화학 조절액 자동 카트리리지, Stat Profile Prime Plus® VET 血气、一氧化二碳-血氧仪、化学对照液自动试剂盒 (含副药)

LOT 24355032

CONTROL 1 2 3 4 5

2026-06-11

Table with columns for Analyte, Unit, and Expected Ranges for Control 1, 2, 3, 4, and 5. Rows include pH, H+, PCO2, PCO2, PO2, PO2, SO2, Hct, Na+, K+, Cl-, iCa, iCa, iMg, Glu, Glu, Lac, BUN, BUN, Urea, Urea, Creatinine, Creatinine, HbF, tHb, tHb, tHb, O2Hb, COHb, MetHb, HHb, tBil, tBil, tBil.

EN

Product Description: Accutest dual-chamber material for monitoring the performance of pH, PCO2, PO2, SO2, hematocrit (Hct), total hemoglobin (tHb), total hemoglobin (tHb), methemoglobin (MetHb), methemoglobin (MetHb), and desoxyhemoglobin (HbF) in levels 1, 2 and 3 as well as Na+, K+, Cl-, Ca, Mg, glucose, lactate, BUN (Urea) and Creatinine in levels 4 and 5. A cartridge of consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation. Part use with Stat Profile Prime Plus VET Analyzers ONLY.

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Intended Use: Intended for monitoring the performance of the Stat Profile Prime Plus VET Analyzers. Controls Levels 1, 2 and 3 are buffered bicarbonate solutions containing O2, salts and preservatives. Each level has a known pH and is equilibrated to a known PCO2, CO2 and N2 value. Controls Level 4 and 5 are buffered solutions containing known concentrations of Na+, K+, Cl-, Ca, Mg, Glucose, Lactate, BUN (Urea), Creatinine and preservatives. Each pouch contains a minimum of 100 mL. Controls contain a known amount of creatinine, however good laboratory practices should be followed during handling of these materials. (REF: NCCLS DOCUMENT M29-72)

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200 Prospect Street Waltham, MA 02454 U.S.A.