REF 57838, 57839, 57840

R_X Only IVD 20°C 8°C CE

Stat Profile Prime Plus® Blood Gas, CO-Oximeter, Chemistry Controls Auto-Cartridge

Cartucho automático para controles de química, gases en sangre y cooxímetro Stat Profile Prime Plus®, Cartouche automatique de contrôles de gaz du sang/CO-oxymètre, chimie Stat Profile Prime Plus®, Stat Profile Prime Plus® Auto-Kassette für Blutgas-, CO-Oximeter- und Blutchemiekontrollen, Stat Profile Prime Plus® Αυτόματο φασίγγιο ελεγχού χημείας αερίων αίματος, CO-Οξυμέτρο, Cartuccia per controlli automatici chimici per gas ematici/CO-ossimetria Stat Profile Prime Plus®, Cartucho automático de controles de química, de CO-oxímetro e de gás no sangue Stat Profile Prime Plus®, Stat Profile Prime Plus® vérgáz/CO-oximéter, Stat Profile Prime Plus® Kan Gaz, CO Oksimetre, Kimya Kontrolleri Otomatik Kartuşu, kémiáitárolókok automatikus patron, Stat Profile Prime Plus® בדיקת גזי דם אוטומטית Co-Oximeter, Stat Profile Prime Plus® 血液ガス, CO オキシメーター, 生化学検査用コントロール自動カートリッジ, Stat Profile Prime Plus® 혈액 가스, CO-산소 농도계, 화학 조절제 자동 카트리지, Stat Profile Prime Plus® 血气、一氧化碳-血氧仪、化学对照溶液自动试剂盒

LOT 24355031

CONTROL 1 2 3 4 5

2026-06-11

EN

Product Description: Aqueous quality control material for monitoring the performance of pH, PCO2, PO2, SO2, hemoglobin (Hct), fetal hemoglobin (HbF), total hemoglobin (HbT), bilirubin total (TBil), urea nitrogen (BUN), creatinine (Cr), aspartate aminotransferase (AST), lactate dehydrogenase (LDH), and aspartate aminotransferase (ALT) in Levels 1, 2 and 3 as well as Na+, K+, Cl-, Ca, Mg, Glucose, and Lactate in Levels 4 and 5. For use with Stat Profile Prime Plus Analyzers ONLY.

Intended Use: Intended for in vitro diagnostic use by healthcare professionals for monitoring the performance of the Stat Profile Prime Plus Analyzers.

Methodology: Refer to Stat Profile Prime Plus Analyzer Instructions For Use Manual for Methodology and Principles.

Composition: Controls Levels 1, 2 and 3 are buffered bicarbonate solutions containing dye, salts and preservatives. Each level has a known pH and is equilibrated to a known O2, CO2, and N2 value. Controls Level 4 and 5 are buffered solutions containing known concentrations of Na+, K+, Cl-, Ca, Mg, Glucose, Lactate and preservatives. Each pouch contains a minimum of 100 mL. Controls contain no constituents of human origin, however good laboratory practices should be followed during handling of these materials. (REF. NCLLS DOCUMENT M29-T2).

Warnings and Cautions: DO NOT FREEZE. Mix the cartridge by gently inverting for several seconds. DO NOT SHAKE CARTRIDGE. Refer to Stat Profile Prime Plus Analyzer Instructions For Use Manual for complete information.

Directions for use: Ensure controls are at room temperature prior to installation. Mix Cartridge well by gently inverting for 1 minute. Verify that the Lot Number on the Expected Ranges Table corresponds to the Lot Number on the cartridge. Refer to Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions.

Limitations: PO2 values vary inversely with temperature (approximately 1%/°C). Therefore, it is critical to follow the temperature guidelines associated in "Directions for Use." The Expected Range values are specific for instruments and controls manufactured by Nova Biomedical. Once installed, each Stat Profile Prime Plus Cartridge may be used for a maximum of 35 days from the initial installation date on the system at which time the system will indicate the cartridge is invalid. Each cartridge may be inserted and removed from the analyzer a maximum of 6 times.

Traceability of Standards: Total Hemoglobin (HbT) and Methemoglobin (MetHb) are traceable by using Cyanmethemoglobin method. Carboxyhemoglobin (COHb) and Oxymethoglobin (OxHb) are traceable using Spectrophotometry. Analytes are traced to NIST Standard Reference Materials.

Reference Intervals: Concentrations are formulated as normal and abnormal expected values in patient blood. The expected clinical range of these values in patient blood is referenced in their, N.W. Hb, 1986 Textbook of Clinical Chemistry, W.B. Saunders Co. Users may wish to determine Mean Values and Expected Ranges in their own laboratory.

Expected Ranges: The expected range for each parameter was determined at Nova Biomedical using replicate determinations on Nova Analyzers. The expected range indicates the maximum deviation from the Mean Value that may be expected under either laboratory conditions for instruments operating within specifications. Refer to Expected Ranges Table.

*Not available in the USA for Point-of-Care/Non-Patient Testing use.

NCLLS Document M29-T2: 2 How to Define and Determine Reference Intervals in the clinical laboratory, approved guidelines-second edition, NCLLS C28-A2, Volume 20, Number 13.

Expected Ranges, Rangos esperados, Plages attendues, Erwartungsbereiche, Αναμενόμενα εύρος, Intervalli previsti, Intervalos previstos, Beklenen Analitlar, Várt tartományok, הצפויים הטווחים, 予測範囲, 예상 범위, 预期范围

Table with 6 columns: Parameter, Unit, CONTROL 1, CONTROL 2, CONTROL 3, CONTROL 4, CONTROL 5. Rows include pH, H+, PCO2, PCO2, PO2, PO2, SO2, Hct, Na+, K+, Cl-, iCa, iCa, iMg, iMg, Glu, Glu, Lac, Lac, HbF*, Hb, Hb, Hb, O2Hb, COHb, MetHb, Hb, TBil*, TBil*, TBil*.

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EC REP Nova Biomedical GmbH Hessenring 13A, Geb. G 64546 Mörfelden-Walldorf Germany
200 Prospect Street Waltham, MA 02454 U.S.A.