
REF 52427,52861,53365,53105,53469,53470,52862,52863,53364,53466,53467,53468

# Stat Profile® Prime Calibrator Cartridge CCS/CCS Comp

Kalibratorkassette, Κασέτα Βαθμονομητή, Cartucho del calibrador, Cartouche étalon, Cartuccia calibratore, Calibrator Cartridge, Kalibratorkassett, Kalibrátorkazetta, キャリブレータカートリッジ, 측정기 카트리지, 校准液盒, Stat Profile® Prime Calibrator Cartridge CCS/CCS Comp

LOT 25015050

2026-06-10

CAL A		≥ 500 mL	CAL B		≥ 250 mL	CAL F		≥ 720 mL	SOLN R		≥ 300 mL
pH	7.345		pH	6.840		PCO <sub>2</sub>	49.5 mmHg		KCl	2 mol/L	
PCO <sub>2</sub>	25.9 mmHg		PO <sub>2</sub>	167.0 mmHg		PO <sub>2</sub>	97.6 mmHg		<b>WASTE CONTAINER</b>	<b>W</b>	
Na <sup>+</sup>	133.0 mmol/L		Na <sup>+</sup>	68.5 mmol/L							
K <sup>+</sup>	4.0 mmol/L		K <sup>+</sup>	10.0 mmol/L							
Cl <sup>-</sup>	106.0 mmol/L		Cl <sup>-</sup>	46.0 mmol/L							
iCa	1.10 mmol/L		iCa	2.20 mmol/L							
Glu	80 mg/dL		Glu	200 mg/dL							
Lac	2.0 mmol/L		Lac	10.0 mmol/L							

**(EN)**  
**Product Description**  
 The above table lists the concentrations of constituents for the calibrators and solutions contained in the Calibrator Cartridge. Contains the following aqueous solutions within individual flexible bags housed in a cardboard box: 3 Calibrators (Calibrator A, B, and F), reference solution, and a flexible waste bag. Each bag includes a filament with a septa. The exposed bag filaments are arranged in a line along the rear of the cartridge box. These septa are pierced during the insertion of the cartridge into the analyzer.

**Intended Use**  
 The Stat Profile Prime Calibrator Cartridge CCS is intended for the calibration of pH, PCO<sub>2</sub>, PO<sub>2</sub>, Hct, Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, iCa, Glucose (Glu), and Lactate (Lac) using the Stat Profile Prime CCS Analyzer.

**Methodology**  
 Refer to Prime Analyzer Instructions For Use Manual for Test Methodology and Principles.

**Composition**  
 Calibrator volumes are listed in the table above. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF: NCCLS DOCUMENT M29-T2).

**Warnings and Cautions:**  
 Intended for in vitro diagnostic use. Refer to the Prime Analyzer Instructions for Use Manual for complete directions for use, installation instructions, recommendations for calibration intervals, troubleshooting information, Methodology and Principles of the testing procedure. Follow standard practices required for handling laboratory reagents. Discard Calibrator Cartridges and Waste products in accordance with local guide lines. Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

**Storage**  
 Store at 15-30°C. Expiration date is printed on the front of the Calibrator Cartridge. Directions For Use Refer to the Prime Analyzer Instructions for Use Manual for complete directions for installation of the cartridge, and recommendations for Calibration Intervals. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use.

**Traceability of Standards**  
 Analytes traced to NIST Standard Reference Materials.

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