
REF 45150, 52864, 53455

Stat Profile® Prime Auto QC Cartridge CCS

Kontroll-Auto-Kassetten, Αυτόματη Κασέτα Διαλυμάτων Ελέγχου, Cartucho automático para controles, Cartouche CQ, Cartuccia automatica controlli, Cartucho automático de control, Kontrollor, autokassett, Automatikus minőség-ellenőrző kazetta, 自動QCカートリッジ, 자동QC카트리지, 自动质控盒, Stat Profile® Prime Auto QC Cartridge CCS

CONTROL 1 2 3

LOT 25010054 2026-06-30

EN Product Description
Consists of 3 flexible bags within a cardboard carton. Each bag contains an aqueous quality control material for monitoring the measurement of pH, PCO₂, PO₂, hematocrit (Hct), Na⁺, K⁺, Cl⁻, Ca, Glucose (Glu) and Lactate (Lac) for use with Nova Biomedical analyzers ONLY.
Formulated at three levels:
CONTROL 1 Acidosis, with High Electrolyte, Low Normal Glu, Normal Lac
CONTROL 2 Normal pH, Low Normal Hct, Normal Electrolyte, High Glu, High Lac
CONTROL 3 Alkalosis, High Hct, Low Electrolyte, High Abnormal Glu, High Abnormal Lac

Intended Use
The StatProfile Prime Auto QC Cartridge CCS is a Quality Control material intended for *in vitro* diagnostic use by healthcare professionals for monitoring the performance of the StatProfile Prime CCS Analyzer.

Methodology
Refer to Prime Analyzer Instructions For Use Manual.

Composition
A buffered bicarbonate solution, each control with a known pH and known levels of Na⁺, K⁺, Cl⁻, Ca, Glu and Lac. Solutions are equilibrated with known levels of O₂, CO₂ and N₂. The conductivity is equivalent to a known hematocrit value in whole blood. Mold inhibited. Each bag contains a minimum volume of 100 mL. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-T2).

Warnings and Cautions:
Must be stored at 24-26°C for at least 24 hours before use. DO NOT FREEZE. Intended for *in vitro* diagnostic use. Refer to the Prime Analyzer Instructions for Use for complete instructions for use, including recommendations for use of controls, troubleshooting information, Methodology and Principles of the testing procedures. Follow standard practices required for handling laboratory reagents.

Storage
Store at 2-8°C. DO NOT FREEZE. Expiration date is printed on each cartridge.

Directions for use
Cartridge must be stored at approximately 24-26°C for at least 24 hours prior to opening. Refer to Analyzer Instructions for Use for complete instructions. Verify that the lot number appearing on the Expected Ranges Table is the same lot number indicated on the control cartridge. The cartridge should be mixed by gently inverting for several seconds. Do not shake the cartridge.

Limitations
pO₂ values vary inversely with temperature (approximately 1%/°C). The Expected Range values are specific for instruments and calibrators manufactured by Nova Biomedical.

Traceability of Standards
Analyses are traced to NIST Standard Reference Materials.

Reference Intervals
Concentrations are formulated to represent three pH levels (Acidosis, Normal pH, and Alkalosis).

The expected clinical range of values for these analytes in patient fluids is defined in the following:
1986 Textbook of Clinical Chemistry, W.B. Saunders Co.
Users may wish to determine MEAN VALUES and EXPECTED RANGE in their own laboratory.

Expected Ranges
The EXPECTED RANGE for each analyte was determined at Nova Biomedical by using multiple runs of each level of control at 37°C on multiple instruments.

The EXPECTED RANGE indicates the maximum deviation from the mean value that may be expected under differing laboratory conditions for instruments operating within specifications.

Refer to Expected Ranges Table.

How to Define and Determine Reference Intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13

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EC REP Nova Biomedical GmbH Hessenring 13A, Geb. G 64546 Mörfelden-Walldorf Germany
Waltham, MA 02454, U.S.A.