Stat Profile Prime Plus® Blood Gas, CO-Oximeter, Chemistry Controls Auto-Cartridge

Cartucho automático para controles de química, gases en sangre y coximetro Stat Profile Prime Plus®, Cartouche automatique de contrôles de gaz du sang/CO-oxymètre, chimie Stat Profile Prime Plus®, Stat Profile Prime Plus® Auto-Kassette für Blutgas-, CO-Oximeter- und Blutchemiekontrollen, Stat Profile Prime Plus® Αυτόματη φασίγγιο εμπίεδο ελέγχου χημείας ερίσμου αίματος, CO-Οξυμέτρο, Cartuccia per controlli automatici chimici per gas ematici/CO-ossimetria Stat Profile Prime Plus®, Cartucho automático de controles de química, de CO-oxímeter e de gás no sangue Stat Profile Prime Plus®, Stat Profile Prime Plus® vėrgėzė/CO-oximėter, Stat Profile Prime Plus® Kan-Gaz, CO Oxsimeter, Kimya Kontrolleri Otomatik Kartuşu, kėmiyal kontrollek otomatikus patron, Stat Profile Prime Plus® מוסמית לוג בדם Co-Oximeter, Stat Profile Prime Plus® 血液ガス, CO オキシメーター, 生化学検査用コントロール自動カートリッジ, Stat Profile Prime Plus® 혈액 가스, CO-산소 농도계, 화학 조절제 자동 카트리지, Stat Profile Prime Plus® 血气、一氧化碳-血氧仪、化学对照溶液自动试剂盒

LOT 25007024

CONTROL 1 2 3 4 5

2026-06-20

EN

Product Description: Accutest quality reagent material for monitoring the performance of pH, PCO2, PO2, SO2, hemato-crit, total hemoglobin (tHb), total hemoglobin (tHb), bicarbonate (HCO3-), lactate, glucose, urea, creatinine, electrolytes (Na+, K+, Cl-, Ca, Mg, glucose, lactate) and osmolarity (Osm) on the Stat Profile Prime Plus Analyzers ONLY.

Intended Use: Intended for use by trained professionals for monitoring the performance of the Stat Profile Prime Plus Analyzers. Refer to the Stat Profile Prime Plus Analyzer Instructions for Use Manual for Verification and Principles.

Composition: Controls 1, 2 and 3 are buffered electrolyte solutions containing electrolyte salts and preservatives. Each level has a known amount of Na+, K+, Cl-, Ca, Mg, glucose, lactate and osmolarity. Each level contains a maximum of 100 mg/L of each analyte. Controls 4 and 5 are unbuffered electrolyte solutions containing electrolyte salts and preservatives. Each level contains a maximum of 100 mg/L of each analyte.

Warnings and Cautions: DO NOT FREEZE. Use the cartridge by gently inverting for several seconds. DO NOT SHAKE CARTRIDGE. Refer to the Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete information.

Storage: Store at 2-8°C (37-46°F). DO NOT FREEZE. Directly for use: Ensure cartridge is at room temperature prior to installation. Vx Cartridge will go gently, inverting for 1 minute. Verify that the Lot Number on the Expected Ranges Table corresponds to the Lot Number on the cartridge. Refer to Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions.

Limitations: PO2 values vary inversely with temperature (approximately 1%/°C). Therefore, it is not safe to follow the temperature guidelines in "Directions for Use." The Expected Range values are specific for instruments and controls manufactured by Nova Biomedical. Once installed, each Stat Profile Prime Plus Cartridge may be used for a maximum of 30 days from the initial installation date on the system (when the system will indicate the cartridge is viable). Each cartridge may be replaced and removed from the analyzer at any time.

Traceability of Standards: Total hemoglobin (tHb) and Methemoglobin (MetHb) are traceable by using Cyanmethemoglobin method. Carboxyhemoglobin (COHb) and Oxyhemoglobin (O2Hb) are traceable using Spectrophotometric analysis as traced to NIST Standard Reference Materials.

Reference Intervals: Concentrations are formulated at normal and abnormal expected values in patient blood. The expected clinical range of these values in patient blood is referenced in Tietz, N.Y. et al. 1986 Textbook of Clinical Chemistry, W.B. Saunders Co. Users may wish to determine their own Reference Intervals and Expected Ranges in their own laboratory.

Expected Ranges: The expected range for each parameter was determined at Nova Biomedical using replicate determinations on Nova analyzers. The expected range indicates the maximum deviation from the Mean Value that may be expected under differing laboratory conditions for instruments operating within specifications. Refer to Expected Ranges Table.

*Not available in the USA for Point-of-Care/Bedside Patient Testing.

INCCLS Document M25-72: Show to Define and Determine Reference Intervals in the clinical laboratory, approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13.

Expected Ranges, Rangos esperados, Plages attendues, Erwartungsbereiche, Αναμενόμενα εύρη, Intervalli previsti, Intervallios previstos, Beklenen Aralıklar, Vårt tartományok, הצפויים הטווחים, 予測範囲, 예상 범위, 预期范围

Table with 6 columns: Parameter, CONTROL 1, CONTROL 2, CONTROL 3, CONTROL 4, CONTROL 5. Rows include pH, H+, PCO2, PO2, SO2, P02, Hct, Na+, K+, Cl-, iCa, iCa, iMg, iMg, Glu, Glu, Lac, Lac, HbF*, tHb, tHb, tHb, O2Hb, COHb, MetHb, HbH, tBil*, tBil*, tBil*.

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