
REF 53107, 53108, 53456

Stat Profile® Prime Auto QC Cartridge ABG

Kontroll-Auto-Kassetten, Αυτόματη Κασέτα Διαλυμάτων Ελέγχου, Cartucho automático para controles, Cartouche CQ, Cartuccia automatica controlli, Cartucho automático de control, Kontroller, autokassett

CONTROL 1 2 3

LOT 25010054 2026-06-30

| | | Expected Ranges, Erwartete Bereiche, Αναμενόμενα Εύρη, Rangos esperados, Fourchettes attendues, Intervalli previsti, Gamas previstas, Förväntade områden | | |
|------------------|--------|--|-----------------------|-----------------------|
| | | CONTROL 1 | CONTROL 2 | CONTROL 3 |
| | | min - \bar{x} - max | min - \bar{x} - max | min - \bar{x} - max |
| pH | | 7.111 - 7.141 - 7.171 | 7.325 - 7.355 - 7.385 | 7.553 - 7.583 - 7.613 |
| H ⁺ | nmol/L | 77.45 - 72.28 - 67.45 | 47.32 - 44.16 - 41.21 | 27.99 - 26.12 - 24.38 |
| pCO ₂ | mmHg | 66.5 - 73.5 - 80.5 | 41.0 - 46.0 - 51.0 | 18.3 - 22.3 - 26.3 |
| pCO ₂ | kPa | 8.87 - 9.80 - 10.73 | 5.47 - 6.13 - 6.80 | 2.44 - 2.97 - 3.51 |
| pO ₂ | mmHg | 49.1 - 59.1 - 69.1 | 90.4 - 100.4 - 110.4 | 130.3 - 145.3 - 160.3 |
| pO ₂ | kPa | 6.55 - 7.88 - 9.21 | 12.05 - 13.39 - 14.72 | 17.37 - 19.37 - 21.37 |

Product Description
Consists of 3 flexible bags within a cardboard carton. Each bag contains an aqueous quality control material for monitoring the measurement of pH, PCO₂, and PO₂ for use with Nova Biomedical analyzers ONLY.

Formulated at three levels:
 Acidosis
 Normal pH
 Alkalosis

Intended Use
The Stat Profile Prime Auto QC Cartridge ABG is a Quality Control Material intended for in vitro diagnostic use by healthcare professionals for monitoring the performance of the Stat Profile Prime ABG Analyzer.

Methodology
Refer to Prime Analyzer Instructions For Use Manual.

Composition
A buffered bicarbonate solution, each control with a known pH. Solutions are equilibrated with known levels of O₂, CO₂, and N₂. Mold inhibited. Each bag contains a minimum volume of 100 mL. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF: NCCLS DOCUMENT M29-T2).

Warnings and Cautions:
Must be stored at 2-26°C for at least 24 hours before use. DO NOT FREEZE. Intended for in vitro diagnostic use. Refer to the Prime Analyzer Instructions for Use for complete directions for use, including recommendations for use of controls, troubleshooting information, Methodology and Principles of the testing procedures. Follow standard practices required for handling laboratory reagents.

Storage
Store at 2-8°C; DO NOT FREEZE. Expiration date is printed on each cartridge.

Directions for use
Cartridge must be stored at approximately 24-26°C for at least 24 hours prior to opening. Refer to Analyzer Instructions for Use for complete instructions. Verify that the lot number appearing on the Expected Ranges Table is the same lot number indicated on the control cartridge. The cartridge should be mixed by gently inverting for several seconds. Do not shake the cartridge.

Limitations
PO₂ values vary inversely with temperature (approximately 1%/°C). The Expected Range values are specific for instruments and calibrators manufactured by Nova Biomedical.

Traceability of Standards
Analyses are traced to NIST Standard Reference Materials.

Reference Intervals
Concentrations are formulated to represent three pH levels (Acidosis, Normal pH, and Alkalosis). The expected clinical range of values for these analytes in patient blood is referenced in Tietz, NW ed 1986 Textbook of Clinical Chemistry, W.B. Saunders Co. Users may wish to determine MEAN VALUES and EXPECTED RANGES in their own laboratory.*

Expected Ranges
The EXPECTED RANGE for each analyte was determined at Nova Biomedical by using multiple runs of each level of control at 37°C on multiple instruments. The EXPECTED RANGE indicates the maximum deviations from the mean value that may be expected under differing laboratory conditions for instruments operating within specifications. Refer to Expected Ranges Table.

*How to Define and Determine Reference Intervals in the clinical laboratory, approved guidelines-second edition, NCCLS C28-A2, Volume 20, Number 13

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