Stat Profile Prime Plus® VET Blood Gas, CO-Oximeter, Chemistry Controls Auto-Cartridge with Creatinine

Cartucho automático con creatinina para controles de química, gases en sangre y cooxímetro Stat Profile Prime Plus® VET, Cartouche automatique de contrôles de gaz du sang/CO-oxymètre, chimie Stat Profile Prime Plus® VET avec créatine, Stat Profile Prime Plus® VET Auto-Kassette mit Kreatinin für Blutgas-, CO-Oximeter- und Blutchemiekontrollen, Αυτόματη συσκευή με κρεατινίνη εμπρόθετος εξοπλισμός χημείας ερίστων αίματος, CO-Οξυμέτρο Stat Profile Prime Plus® VET, Cartucucia con creatinina per controlli automatici chimici per gas ematici/CO-ossimetria Stat Profile Prime Plus® VET, Cartucho automático de controles de química, de CO-oxímetro e de gás no sangue Stat Profile Prime Plus® VET com creatinina, Stat Profile Prime Plus® VET vérőgáz, CO-oximéter, kémiai kontrollok automatikus patron kreatininnal, (מחסנית בקרת כימיה אוטומטית עם קריאטינין) Stat Profile Prime Plus® VET Chemistry Controls Auto-Cartridge, Stat Profile Prime Plus® VET 血液ガス, CO オキシメーター, 生化学検査用コントロール自動カートリッジ (クレアチニン), 크레아티닌 사용 Stat Profile Prime Plus® VET 혈액 가스, CO-산소 농도계, 화학 조절제 자동 카트리지, Stat Profile Prime Plus® VET 血气、一氧化炭-血氧仪、化学对照液自动试剂盒 (含肌酐)

LOT 25007025

2026-06-20

CONTROL 1 2 3 4 5

Table with 6 columns: Parameter, Unit, CONTROL 1 (min - X - max), CONTROL 2 (min - X - max), CONTROL 3 (min - X - max), CONTROL 4 (min - X - max), CONTROL 5 (min - X - max). Rows include pH, H+, Pco2, Pcoo2, PO2, PO2, SO2, Hct, Na+, K+, Cl-, iCa, img, img, Glu, Glu, Lac, Lac, BUN, BUN, Urea, Urea, Creatinine, Creatinine, Creatinine, HbF, iHb, iHb, O2Hb, COHb, MetHb, HbH, iBil, iBil, iBil.

EN

Product Description: Accurate quality control material for monitoring the performance of pH, Pco2, Pcoo2, SO2, hematocrit (Hct), total hemoglobin (tHb), total hemoglobin (tHb), hemoglobin (Hb), carboxyhemoglobin (COHb), and desoxyhemoglobin (dHb) in Levels 1, 2 and 3 as well as Na+, K+, Cl-, Ca, Mg, Glucose, Lactate, BUN (urea), and Creatinine in Levels 4 and 5. A charge kit consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation. For use with Stat Profile Prime Plus VET Analyzers ONLY.

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Directions for use: Ensure controls are room temperature prior to installation. Level 4 and Level 5 must be charged with creatinine syringes prior to installation of the Calibrator Cartridge. The 2 syringes are labeled and color coded to correspond to the fittings on the cartridge. Activate the cartridge as follows: 1. Hold the syringe with top side down and remove protective cap. 2. Attach one of the enclosed needle assemblies to the syringe. Remove the protective cover from the needle. 3. Match the color and label of the syringe to the appropriate fitting and insert needle. 4. Slowly depress syringe plunger until the contents are dispensed. DO NOT PULL BACK ON THE PLUNGER TO FLUSH CONTENTS OF SYRINGE. 5. Repeat steps 1-4 for the next Control. 6. Repeat Steps 1-4 for the next Control. 7. Make sure the cartridge is completely covered within 1 minute. Cartridge is ready for use. Verify that the Lot Number on the Expected Ranges Table corresponds to the Lot Number on the cartridge. Refer to Stat Profile Prime Plus VET Analyzer Instructions for Use Manual for complete directions.

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200 Prospect Street Waltham, MA 02454 U.S.A.