
REF 57831, 57832, 57833, 57834, 57835, 65142

Stat Profile Prime Plus® Calibrator Cartridge With Creatinine

Stat Profile Prime Plus® Kalibrierkassette mit Kreatinin, Stat Profile Prime Plus® Φυσίγγιο μέσου βαθμονόμησης με κρεατινίνη, Cartucho de calibrador con creatinina Stat Profile Prime Plus®, Cartouche de calibreurs Stat Profile Prime Plus® avec créatinine, Cartuccia calibratore con creatinina Stat Profile Prime Plus®, Cartucho de calibrador Stat Profile Prime Plus® com creatinina, Stat Profile Prime Plus® kalibrátor patron kreatininne, מוסכיתת הכויל עם קריאטינין, Stat Profile Prime Plus®, Stat Profile Prime Plus® キャリブレータカートリッジ (クレアチニン), 크레아티닌 사용 Stat Profile Prime Plus® 보정기 카트리지, Stat Profile Prime Plus® 含肌酐的校准试剂盒, Stat Profile Prime Plus® Kreatinin İçeren Calibrator Cartridge

LOT 25030034

2026-07-21

CAL A		CAL B		CAL C		CAL D		CAL F		SOLN R	
HCO ₃	11 mmol/L	pH	6.834 @ 37°C	HCO ₃	22 mmol/L	pH	7.354 @ 37°C	PO ₂	103.0 mmHg	KCl	2 mol/L
Na ⁺	135.5 mmol/L	H ⁺	146.55 nmol/L	Cl ⁻	69.0 mmol/L	H ⁺	44.26 nmol/L	PO ₂	13.70 kPa	 WASTE CONTAINER	
K ⁺	4.0 mmol/L	Na ⁺	77.0 mmol/L	Glu	200 mg/dL	Glu	11.10 mmol/L				
Cl ⁻	108.0 mmol/L	K ⁺	10.0 mmol/L	Glu	11.10 mmol/L	Lac	10.0 mmol/L				
iCa	1.04 mmol/L	iCa	1.92 mmol/L	Lac	10.0 mmol/L	Lac	89.0 mg/dL				
iCa	4.16 mg/dL	iCa	7.68 mg/dL	BUN	50 mg/dL	BUN	17.9 mmol/L				
iMg	0.50 mmol/L	iMg	1.50 mmol/L	BUN	50 mg/dL	Creat	5.0 mg/dL				
iMg	1.22 mg/dL	iMg	3.63 mg/dL	Creat	5.0 mg/dL	Creat	0.44 mmol/L				
Glu	80 mg/dL	tHb	15.4 g/dL	Creat	0.44 mmol/L	Creat	442 µmol/L				
Glu	4.44 mmol/L	tHb	154 g/L								
Lac	2.0 mmol/L	tHb	9.6 mmol/L								
Lac	17.8 mg/dL										
BUN	10 mg/dL										
BUN	3.6 mmol/L										
Creat	1.0 mg/dL										
Creat	0.09 mmol/L										
Creat	88.4 µmol/L										

EN Description
The Cartridge contains the following aqueous solutions within individual bags: 5 Calibrators (Calibrator A, B, C, D and F), Reference Solution, and a flexible waste bag. Each bag includes a filament with a septum. The exposed bag filaments are arranged in a line along the rear of the cartridge box. The septa are pierced during the insertion of the cartridge into the analyzer. A charge kit consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation.

Intended Use
For the quantitative determination of Na⁺, K⁺, Cl⁻, iCa, iMg, BUN (Urea), Glucose, Lactate, Creatinine, pH, PCO₂, PO₂, SO₂%, Hct, Hb, HbF, O₂Hb, COHb, MetHb, HbB, and tHb in human blood using the Stat Profile Prime Plus Analyzer.

Methodology
Refer to Stat Profile Prime Plus Analyzer Instructions For Use Manual for Methodologies and Principles.

Composition
The calibrators are aqueous salt solutions containing analytes at concentrations listed in the table above. Calibrator B also contains dye. Contain no constituents of human origin, however, good laboratory practice should be followed during the handling of these materials. (REF: NCCLS DOCUMENT M29-T2).

Warnings and Cautions
Intended for *in vitro* diagnostic use.
Refer to the Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions. Follow standard procedure for handling laboratory reagents. Discard Cartridge in accordance with local guidelines.
Visable Bag Contains Preservative (eye, skin, and respiratory irritant).

Storage
Store at 15 - 30°C. Expiration date is printed on the Calibrator Cartridge.

Directions For Use
Cal A and Cal C must be charged with creatinine syringe prior to installation of the Calibrator Cartridge. The two syringes are labeled and color coded to correspond to the filaments on the carton for Cal A and Cal C.

- Hold the syringe with tip side down and remove protective cap.
- Attach one of the enclosed needle assemblies to the syringe. Remove the protective cover from the needle.
- Match the color and label of the syringe to the appropriate filament and insert needle.
- Slowly depress syringe plunger until the contents are dispensed.

DO NOT PULL BACK ON THE PLUNGER TO FLUSH CONTENTS OF SYRINGE.

- Remove needle/syringe assembly from filament and discard in an appropriate sharps container.
- Repeat Steps 1-5 for the next Cal.
- Mix Cartridge well by gently shaking for 1 minute.

Refer to the Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use.

Not available in the USA for Point-of-Care/Near-Patient Testing use.

Traceability of Standards
Analytes traced to NIST Standard Reference Materials wherever possible.

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EC REP Nova Biomedical GmbH Hessenring 13A, Geb. G 64546 Mörfelden-Walldorf Germany
Waltham, MA 02454 U.S.A.