REF 48903, 48904

# Stat Profile® pH/Ox Ultra/Critical Care Xpress Blood Gas Controls Auto-Cartridge

Blutgas-Kontroll-Auto-Kassetten, Αυτόματη Κασέτα Διαλυμάτων Ελέγχου Αερίων Αίματος, Cartucho automático para controles de gases en sangre, Cartouche CQ pour les gaz du sang, Cartuccia automatica controlli, Cartucho automático de controle de gases no sangue, Blodgaskontroller, autokasset

CONTROL 1 2 3

LOT 25034072 2026-07-28

	Expected Ranges, Erwartete Bereiche, Αναμενόμενα Εύρη Rangos esperados, Fourchettes attendues, Intervalli previsti, Gamas previstas, Förväntade områden		
	CONTROL 1	CONTROL 2	CONTROL 3
	min - $\bar{x}$ - max	min - $\bar{x}$ - max	min - $\bar{x}$ - max
pH	7.168 - 7.193 - 7.218	7.369 - 7.394 - 7.419	7.522 - 7.549 - 7.576
H <sup>+</sup>	nmol/L 67.99 - 64.19 - 60.60	42.75 - 40.35 - 38.10	30.09 - 28.28 - 26.57
pCO <sub>2</sub>	mmHg 54.9 - 60.2 - 65.4	40.7 - 44.2 - 47.7	21.3 - 23.8 - 26.8
pCO <sub>2</sub>	kPa 7.30 - 8.00 - 8.69	5.41 - 5.88 - 6.34	2.84 - 3.17 - 3.57
pO <sub>2</sub>	mmHg 44.9 - 52.9 - 60.9	97.1 - 105.1 - 113.1	123.3 - 135.3 - 147.3
pO <sub>2</sub>	kPa 5.97 - 7.03 - 8.10	12.92 - 13.98 - 15.05	16.40 - 17.99 - 19.59
SO <sub>2</sub>	% 54.6 - 59.6 - 63.6		88.7 - 92.7 - 96.7
Hct	% 32 - 34 - 37		46 - 49 - 52
Hb	g/dL 9.6 - 11.1 - 12.6		14.8 - 16.3 - 17.8
Hb	mmol/L 5.95 - 6.89 - 7.82		9.17 - 10.10 - 11.03

**Product Description**  
An aqueous quality control material for monitoring the measurement of pH, pCO<sub>2</sub>, pO<sub>2</sub>, SO<sub>2</sub>, hematoct (Hct), and hemoglobin (Hb) for use with Nova Biomedical analyzers ONLY. Each cartridge contains controls formulated at three clinically significant levels.

- CONTROL 1 Acidosis, with low SO<sub>2</sub>, low-normal Hct/Hb
- CONTROL 2 Normal pH
- CONTROL 3 Alkalosis, with SO<sub>2</sub>, High Hct/Hb

**Intended Use**  
For in vitro diagnostic use for monitoring the performance of Nova Biomedical Stat Profile pH/Ox Ultra/Critical Care Xpress Analyzers.

**Methodology**  
Refer to Stat Profile pH/Ox Ultra/CCX Analyzer Instructions For Use Manual.

**Composition**  
A buffered bicarbonate solution, the Controls Auto-Cartridge consists of 3 flexible bags, each with a known pH. Solutions are equilibrated with known levels of O<sub>2</sub>, CO<sub>2</sub>, and H<sub>2</sub>. The reflective characteristics are equivalent to a known hemoglobin value in whole blood. The reflective and conductivity signal is equivalent to a known hemoglobin value in blood. (Mold indicated. Each control flexible pouch contains a minimum value of 100mL. Contains no constituents of human origin, however good laboratory practice should be followed during handling of these materials. (REF: NCCLS LAC/UR/INT M29-T2).

**Warnings and Cautions:**  
Must be stored at 2-26°C for at least 24 hours before using. DO NOT SHAKE the cartridge. Intended for in vitro diagnostic use. Refer to the Stat Profile pH/Ox Ultra/Critical Care Xpress Analyzer Instructions for Use for complete directions for use, recommendations for use of controls, troubleshooting information, and methodology and Principles of the testing procedures. Follow standard practices required for handling laboratory reagents.

**Storage**  
Store at 2-26°C. DO NOT FREEZE. It is critical to follow the room temperature equilibrium instructions prior to use as described in "Directions for Use". Each cartridge has a lot number and expiration date printed on the label.

**Directions for Use**  
Controls must be stored at approximately 2-26°C for at least 24 hours prior to opening. Refer to Analyzer Instructions for Use for complete instructions. Verify that the lot number appearing on the Expected Ranges Table is the same lot number indicated on the control material being analyzed. The Auto-Cartridge should be mixed by gently inverting the pack for several seconds. DO NOT SHAKE the cartridge.

**Limitations**  
pO<sub>2</sub> values vary inversely with temperature (approximately 1%/°C). The Expected Range values are specific for instruments and controls manufactured by Nova Biomedical. Once installed, each Stat Profile pH/Ox Ultra/CCX Auto-Cartridge may be used for a maximum of 35 days from the initial installation date on the system at which time the system will indicate the cartridge is invalid.

Each Stat Profile pH/Ox Ultra/Critical Care Xpress Auto-Cartridge may be inserted and removed from the pH/Ox Ultra/CCX analyzer a maximum of six times.

**Traceability of Standards**  
Analyses are traced to NIST Standard Reference Materials.

**Reference Intervals**  
Concentrations are formulated to represent three pH levels (Acidosis, Normal pH, and Alkalosis), and with Low-Normal and High Hct/Hb levels.

The expected clinical range of values for these analytes in patient blood is referenced in Tietz, NW ed 1996 Textbook of Clinical Chemistry, W.B. Saunders Co.

Users may wish to determine MEAN VALUES and EXPECTED RANGES in their own laboratory!

**Expected Ranges**  
The EXPECTED RANGE for each analyte was determined at Nova Biomedical by using multiple runs of each level of control at 37°C on multiple instruments. The EXPECTED RANGE indicates the maximum deviation from the mean value that may be expected under differing laboratory conditions for instruments operating within specifications. Refer to Expected Ranges Table.

\*How to Define and Determine Reference Intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13

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