Stat Profile® Prime VET CCS Comp Calibrator Cartridge

Kalibratorkassette, Κάσέτα Βαθμονομητή, Cartucho del calibrador, Cartouche étalon, Cartuccia calibratore, Calibrator Cartridge, Kalibratorkassett, Kalibrátorkazetta, 교정기 카트리지, キャリブレータカートリッジ, 校准液盒

LOT 25065054

2026-08-24

CAL A ≥ 500 mL		CAL B ≥ 250 mL		CAL F ≥ 720 mL		SOLN R ≥ 300 mL	
pH	7.354	pH	6.840	PCO ₂	50.7 mmHg	KCl	2 mol/L
PCO ₂	24.7 mmHg	PO ₂	167.0 mmHg	PO ₂	101.1 mmHg	 WASTE CONTAINER W	
Na ⁺	133.0 mmol/L	Na ⁺	68.5 mmol/L				
K ⁺	4.0 mmol/L	K ⁺	10.0 mmol/L				
Cl ⁻	106.0 mmol/L	Cl ⁻	46.0 mmol/L				
iCa	1.10 mmol/L	iCa	2.20 mmol/L				
Glu	80 mg/dL	Glu	200 mg/dL				
Lac	2.0 mmol/L	Lac	10.0 mmol/L				

EN

Product Description
The above table lists the concentrations of constituents for the calibrators and solutions contained in the Calibrator Cartridge. Contains the following aqueous solutions within individual flexible bags housed in a cardboard box: 3 Calibrators (Calibrator A, B, and F), reference solution, and a flexible waste bag. Each bag includes a firmet with a septa. The exposed bag firmets are arranged in a line along the rear of the cartridge box. These septa are pierced during the insertion of the cartridge into the analyzer.

Intended Use
For the quantitative determination of Na, K, Cl, iCa, Glucose, Lactate, pH, PCO₂, PO₂, and Hct in blood using the Stat Profile Prime VET CCS Comp Analyzer.

Methodology
Refer to Prime Vet CCS Comp Analyzer Instructions For Use Manual for Test Methodologies and Principles.

Composition
CALibrator volumes are listed in the table above. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-T2).

Warnings and Cautions:
For Veterinary Use Only, Not for Human Use. Refer to the Prime Vet CCS Comp Analyzer Instructions for Use Manual for complete directions for use, installation instructions, and recommendations for calibration intervals, troubleshooting information, Methodology and Principles of the testing procedures. Follow standard practices required for handling laboratory reagents. Discard Calibrator Cartridges and Waste products in accordance with local guide lines. Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

Storage
Store at 15-30°C. Expiration date is printed on the front of the Calibrator Cartridge.

Directions For Use
Refer to the Prime Vet CCS Comp Analyzer Instructions for Use Manual for complete directions for installation of the cartridge, and recommendations for Calibration Intervals. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use.

Traceability of Standards
Analytes traced to NIST Standard Reference Materials.

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