# Stat Profile Prime Plus<sup>®</sup> VET Calibrator Cartridge With Creatinine

Stat Profile Prime Plus<sup>®</sup> VET Kalibriermessung mit Kreatinin, Φασίγγιο μέσου βαθμονόμησης VET με κρεατινίνη Stat Profile Prime Plus<sup>®</sup>, Cartucho de calibrador con creatinina Stat Profile Prime Plus<sup>®</sup> VET, Cartouche de calibreuses avec creatinine Stat Profile Prime Plus<sup>®</sup> VET, Cartuccia calibratore con creatinina Stat Profile Prime Plus<sup>®</sup> VET, Cartucho de calibrador Stat Profile Prime Plus<sup>®</sup> VET con creatinina, Stat Profile Prime Plus<sup>®</sup> VET kalibrátor patron kreatininnel, Stat Profile Prime Plus<sup>®</sup> VET, Cartucho de calibrador Stat Profile Prime Plus<sup>®</sup> VET con creatinina, Stat Profile Prime Plus<sup>®</sup> VET kalibrátor patron kreatininnel, Stat Profile Prime Plus<sup>®</sup> VET קריאטינין כיוול מוחסנית, Stat Profile Prime Plus<sup>®</sup> VET キャリブレタ(クレアチニン), 크레아티닌 사용 Stat Profile Prime Plus<sup>®</sup> VET 보정기 카트리지, Stat Profile Prime Plus<sup>®</sup> VET 含肌酐的校准试剂盒

LOT 25069072

2026-08-12

CAL A ≥ 495 mL		CAL B ≥ 250 mL		CAL C ≥ 245 mL		CAL D ≥ 250 mL		CAL F ≥ 850 mL		SOLN R ≥ 500 mL	
HCO <sub>3</sub>	11 mmol/L	pH	6.834 @ 37°C	HCO <sub>3</sub>	22 mmol/L	pH	7.354 @ 37°C	PO <sub>2</sub>	103.8 mmHg	KCl	2 mol/L
Na <sup>+</sup>	135.5 mmol/L	H <sup>+</sup>	146.55 nmol/L	Cl <sup>-</sup>	69.0 mmol/L	H <sup>+</sup>	44.26 nmol/L	PO <sub>2</sub>	13.81 kPa	<div style="border: 1px solid black; padding: 5px; text-align: center;"> <b>WASTE CONTAINER</b> </div> <div style="border: 1px solid black; padding: 5px; text-align: center; margin-top: 5px;"> <b>W</b> </div>	
K <sup>+</sup>	4.0 mmol/L	Na <sup>+</sup>	77.0 mmol/L	Glu	200 mg/dL						
Cl <sup>-</sup>	108.0 mmol/L	K <sup>+</sup>	10.0 mmol/L	Glu	11.10 mmol/L						
iCa	1.04 mmol/L	iCa	1.92 mmol/L	Lac	10.0 mmol/L						
iCa	4.16 mg/dL	iCa	7.68 mg/dL	Lac	89.0 mg/dL						
iMg	0.50 mmol/L	iMg	1.50 mmol/L	BUN	50 mg/dL						
iMg	1.22 mg/dL	iMg	3.63 mg/dL	BUN	17.9 mmol/L						
Glu	80 mg/dL	tHb	15.6 g/dL	Creat	5.0 mg/dL						
Glu	4.44 mmol/L	tHb	156 g/L	Creat	0.44 mmol/L						
Lac	2.0 mmol/L	tHb	9.7 mmol/L	Creat	442 µmol/L						
Lac	17.8 mg/dL										
BUN	10 mg/dL										
BUN	3.6 mmol/L										
Creat	1.0 mg/dL										
Creat	0.09 mmol/L										
Creat	88.4 µmol/L										

**EN Description**  
The Cartridge contains the following aqueous solutions within individual bags: 5 Calibrators (Calibrator A, B, C, D and F), Reference Solution, and a flexible waste bag. Each bag includes a filament with a septum. The exposed bag filaments are arranged in a line along the rear of the cartridge box. The septa are pierced during the insertion of the cartridge into the analyzer. A charge kit consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation.

**Intended Use**  
For the quantitative determination of Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, iCa, iMg, BUN (Urea), Glucose, Lactate, Creatinine, pH, PCO<sub>2</sub>, PO<sub>2</sub>, SO<sub>2</sub>%, Hct, Hb, O<sub>2</sub>Hb, COHb, MeHb, Hb, and tHb in animal whole blood using the Stat Profile Prime Plus VET Analyzer.

**Methodology**  
Refer to Stat Profile Prime Plus VET Analyzer Instructions For Use Manual for Methodologies and Principles.

**Composition**  
Calibrator and solution volumes are listed in the table. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-T2).

**Warnings and Cautions**  
Refer to the Stat Profile Prime Plus VET Analyzer Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory reagents. Discard Cartridges in accordance with local guidelines. Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

**Storage**  
Store at 15 - 30°C. Expiration date is printed on the Calibrator Cartridge.

**Directions For Use**  
**Activation of Cartridge**  
Cal A and Cal C must be charged with creatinine syringe prior to installation of the Calibrator Cartridge. The two syringes are labeled and color coded to correspond to the filaments on the carton for Cal A and Cal C.  
 1. Hold the syringe with tip side down and remove protective cap.  
 2. Attach one of the enclosed needle assemblies to the syringe. Remove the protective cover from the needle.  
 3. Match the color and label of the syringe to the appropriate filament and insert needle.  
 4. Slowly depress syringe plunger until the contents are dispensed.  
 DO NOT PULL BACK ON THE PLUNGER TO FLUSH CONTENTS OF SYRINGE.  
 5. Remove needle/syringe assembly from filament and discard in an appropriate sharps container.  
 6. Repeat Steps 1-5 for the next Cal.  
 7. Mix Cartridge well by gently shaking for 1 minute.

Refer to the Stat Profile Prime Plus VET Analyzer Instructions for Use Manual for complete directions. Calibrator Cartridge must be inserted gently several times to thoroughly mix contents. Cartridge is ready for use.

**Traceability of Standards**  
Analyses traced to NIST Standard Reference Materials wherever possible.

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