Stat Profile® Prime Calibrator Cartridge CCS/CCS Comp

Kalibratorkassette, Κασέτα Βαθμονομητή, Cartucho del calibrador, Cartouche étalon, Cartuccia calibratore, Calibrator Cartridge, Kalibratorkassett, Kalibratorkazetta, キャリブレータ カートリッジ, 측정기 카트리지, 校准液盒, Stat Profile® Prime Calibrator Cartridge CCS/CCS Comp

LOT 25073039 2026-09-03

| CAL A ≥ 500 mL | | CAL B ≥ 250 mL | | CAL F ≥ 720 mL | | SOLN R ≥ 300 mL | |
|------------------|--------------|-----------------|-------------|------------------|------------|-----------------------------------------------------------------------------------------------------------------|---------|
| pH | 7.355 | pH | 6.840 | PCO ₂ | 49.6 mmHg | KCl | 2 mol/L |
| PCO ₂ | 24.9 mmHg | PO ₂ | 167.0 mmHg | PO ₂ | 100.4 mmHg |  WASTE CONTAINER W | |
| Na ⁺ | 133.0 mmol/L | Na ⁺ | 68.5 mmol/L | | | | |
| K ⁺ | 4.0 mmol/L | K ⁺ | 10.0 mmol/L | | | | |
| Cl ⁻ | 106.0 mmol/L | Cl ⁻ | 46.0 mmol/L | | | | |
| iCa | 1.10 mmol/L | iCa | 2.20 mmol/L | | | | |
| Glu | 80 mg/dL | Glu | 200 mg/dL | | | | |
| Lac | 2.0 mmol/L | Lac | 10.0 mmol/L | | | | |

EN
Product Description
 The above table lists the concentrations of constituents for the calibrators and solutions contained in the Calibrator Cartridge. Contains the following aqueous solutions within individual flexible bags housed in a cardboard box: 3 Calibrators (Calibrator A, B, and F), reference solution, and a flexible waste bag. Each bag includes a ferment with a septa. The exposed bag fittings are arranged in a line along the rear of the cartridge box. These septa are pierced during the insertion of the cartridge into the analyzer.

Intended Use
 The Stat Profile Prime Calibrator Cartridge CCS is intended for the calibration of pH, PCO₂, PO₂, Hct, Na⁺, K⁺, Cl⁻, Glucose (Glu), and Lactate (Lac) using the Stat Profile Prime CCS Analyzer.

Methodology
 Refer to Prime Analyzer Instructions For Use Manual for Test Methodologies and Principles.

Composition
 Calibrator volumes are listed in the table above. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-T2).

Warnings and Cautions:
 Intended for in vitro diagnostic use.
 Refer to the Prime Analyzer Instructions for Use Manual for complete directions for use, installation instructions, recommendations for calibration intervals, Troubleshooting information, Methodology and Principles of the testing procedures. Follow standard practices required for handling laboratory reagents. Discard Calibrator Cartridges and Waste products in accordance with local guide lines. Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

Storage
 Store at 15-30°C. Expiration date is printed on the front of the Calibrator Cartridge. Directions For Use Refer to the Prime Analyzer Instructions for Use Manual for complete directions for installation of the cartridge, and recommendations for Calibration Intervals. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use.

Traceability of Standards
 Analytes traced to NIST Standard Reference Materials.

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