REF 48903, 48904

IVD 20°C $\pm$ 8°C CE

Stat Profile® pHox Ultra/Critical Care Xpress Blood Gas Controls Auto-Cartridge

Blutgas-Kontroll-Auto-Kassetten, Αυτόματη Κασέτα Διαλυμάτων Ελέγχου Αερίων Αίματος, Cartucho automático para controles de gases en sangre, Cartouche CQ pour les gaz du sang, Cartuccia automatica controllati, Cartucho automático de controle de gases no sangue, Blodgaskontrollert, autokassett

CONTROL 1 2 3

LOT 25071033 2026-08-28

Expected Ranges, Erwartete Bereiche, Αναμενόμενα Εύρη Rangos esperados, Fourchettes attendues, Intervalli previsti, Gamas previstas, Förväntade områden

		CONTROL 1	CONTROL 2	CONTROL 3
		min - $\bar{x}$ - max	min - $\bar{x}$ - max	min - $\bar{x}$ - max
pH		7.162 - 7.187 - 7.212	7.375 - 7.400 - 7.425	7.516 - 7.543 - 7.570
H ⁺	nmol/L	68.87 - 65.02 - 61.38	42.21 - 39.85 - 37.62	30.50 - 28.67 - 26.94
pCO ₂	mmHg	57.0 - 62.3 - 67.5	39.9 - 43.4 - 46.9	21.5 - 24.0 - 27.0
pCO ₂	kPa	7.58 - 8.28 - 8.97	5.31 - 5.78 - 6.24	2.85 - 3.19 - 3.59
pO ₂	mmHg	46.7 - 54.7 - 62.7	97.4 - 105.4 - 113.4	126.5 - 138.5 - 150.5
pO ₂	kPa	6.21 - 7.28 - 8.34	12.95 - 14.01 - 15.08	16.83 - 18.42 - 20.02
SO ₂	%	54.2 - 59.2 - 63.2		89.1 - 93.1 - 97.1
Hct	%	32 - 34 - 37		46 - 49 - 52
Hb	g/dL	9.6 - 11.1 - 12.6		14.8 - 16.3 - 17.8
Hb	mmol/L	5.94 - 6.88 - 7.81		9.20 - 10.14 - 11.07

Product Description

An aqueous quality control material for monitoring the measurement of pH, pCO₂, pO₂, SO₂, hemoglobin (Hb), and hemoglobin (Hb) for use with Nova Biomedical analyzers ONLY. Each cartridge contains controls formulated at three clinically significant levels.

CONTROL 1 Acidosis, with low SO₂, low-normal Hct/Hb
CONTROL 2 Normal pH-Wert
CONTROL 3 Alkalosis, with SO₂, High Hct/Hb

Intended Use

For in vitro diagnostic use for monitoring the performance of Nova Biomedical Stat Profile pHox Ultra/Critical Care Xpress Analyzers.

Methodology

Refer to Stat Profile pHox Ultra/CCX Analyzer Instructions for Use Manual.

Composition

A buffered bicarbonate solution, the Controls Auto-Cartridge consists of 3 flexible bags, each with a known pH. Solutions are equilibrated with known levels of O₂, CO₂, and N₂. The resistance characteristics give a signal that is equivalent to a known oxygen saturation value in whole blood. The conductivity characteristics are equivalent to a known hemoglobin value in whole blood. The reflectance and conductivity signal is equivalent to a known hemoglobin value in blood. Mold inhibited. Each control flexible pouch contains a minimum value of 100mls. Contains no constituents of human origin, however good laboratory practice should be followed during handling of these materials. (REF. NCCLS DOCUMENT M23-T2)

Warnings and Cautions:

Must be stored at 24-26°C for at least 24 hours before using. DO NOT SHAKE the cartridge. Intended for in vitro diagnostic use. Refer to the Stat Profile pHox Ultra/CCX Analyzer Instructions for Use for complete directions for use, recommendations for use of controls, troubleshooting information, and methodology and Principles of the testing procedures. Follow standard practices required for handling laboratory reagents.

Storage

Store at 2-8°C. DO NOT FREEZE. It is critical to follow the room temperature equilibrium instructions prior to use as described in "Directions for Use". Each cartridge has a lot number and expiration date printed on the label.

Directions for use

Controls must be stored at approximately 24-26°C for at least 24 hours prior to opening. Refer to analyzer instructions for use for complete instructions. Verify that the lot number appearing on the expected Ranges Table is the same lot number indicated on the control material being analyzed. The Auto-Cartridge should be mixed by gently inverting the pack for several seconds. DO NOT SHAKE the cartridge.

Intimations

Q₁₀ values vary inversely with temperature (approximately 1%/°C).

The Expected Range values are specific for instruments and controls manufactured by Nova Biomedical. Once installed, each Stat Profile pHox Ultra/CCX Auto-Cartridge may be used for a maximum of 35 days from the initial installation date on the system at which time the system will indicate the cartridge is invalid.

Each Stat Profile pHox Ultra/Critical Care Xpress Auto-Cartridge may be inserted and removed from the pHox Ultra/CCX analyzer a maximum of six times.

Traceability of Standards

analyses are traced to NIST Standard Reference Materials.

Reference Intervals

Concentrations are formulated to represent three pH levels (Acidosis, Normal pH, and Alkalosis), and with Low-Normal and High Hct/Hb levels.

The expected clinical range of values for these analytes in patient blood is referenced in Tietz, NW ed 986 Textbook of Clinical Chemistry, W.B. Saunders Co.

Users may wish to determine MEAN VALUES and EXPECTED RANGES in their own laboratory.

Expected Ranges

The EXPECTED RANGE for each analyte was determined at Nova Biomedical using multiple runs of each level of control at 37°C on multiple instruments. The EXPECTED RANGE indicates the maximum deviations from the mean value that may be expected under differing laboratory conditions for instruments operating within specifications. Refer to Expected Ranges Table.

How to Define and Determine Reference Intervals in the clinical laboratory: approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13

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