Stat Profile® Prime Calibrator Cartridge ABG

Kalibratorkassette, Κασέτα Βαθμονομητή, Cartucho del calibrador, Cartouche étalon, Cartuccia calibratore, Calibrator Cartridge, Kalibratorkassett

LOT 25105033

2026-09-18

CAL A ≥ 500 mL pH 7.354 PCO₂ 25.0 mmHg	CAL B ≥ 250 mL pH 6.840 PO₂ 167.0 mmHg	CAL F ≥ 720 mL PCO₂ 50.7 mmHg PO₂ 99.8 mmHg	SOLN R ≥ 300 mL KCl 2 mol/L WASTE CONTAINER W
---	---	---	--

Product Description
The above table lists the concentrations of constituents for the calibrators and solutions contained in the Calibrator Cartridge. Contains the following aqueous solutions within individual flexible bags housed in a cardboard box: 3 Calibrators (Calibrator A, B, and F), reference solution, and a flexible waste bag. Each bag includes a fitment with a septa. The exposed bag fitments are arranged in a line along the rear of the cartridge box. These septa are pierced during the insertion of the cartridge into the analyzer.

Intended Use
For the quantitative determination of pH, PCO₂, and PO₂ in human blood using the Nova Biomedical Prime ABG Analyzer.

Methodology
Refer to Prime Analyzer Instructions For Use Manual for Test Methodologies and Principles.

Composition
Calibrator volumes are listed in the table above. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF: NCCLS DOCUMENT M29-T2).

Warnings and Cautions:
Intended for *in vitro* diagnostic use.
Refer to the Prime Analyzer Instructions for Use Manual for complete directions for use, installation instructions, recommendations for calibration intervals, Troubleshooting information, Methodology and Principles of the testing procedures. Follow standard practices required for handling laboratory reagents. Discard Calibrator Cartridges and Waste products in accordance with local guide lines. Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

Storage
Store at 15-30°C. Expiration date is printed on the front of the Calibrator Cartridge.
Directions For Use Refer to the Prime Analyzer Instructions for Use for complete directions for installation of the cartridge, and recommendations for Calibration Intervals. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use.

Traceability of Standards
Analytes traced to NIST Standard Reference Materials.

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