Stat Profile Prime Plus® VET Blood Gas, CO-Oximeter, Chemistry Controls Auto-Cartridge with Creatinine

Cartucho automático con creatinina para controles de química, gases en sangre y cooxímetro Stat Profile Prime Plus® VET, Cartouche automatique de contrôles de gaz du sang/CO-oxymètre, chimie Stat Profile Prime Plus® VET avec créatine, Stat Profile Prime Plus® VET Auto-Kassette mit Kreatinin für Blutgas-, CO-Oximeter- und Blutchemiekontrollen, Αυτόματο φουίγγιο με κρεατινίνη επιπέδου ελέγχου χημείας αερίων αίματος, CO-Οξυμέτρο Stat Profile Prime Plus® VET, Cartuccia con creatinina per controlli automatici chimici per gas ematici/CO-ossimetria Stat Profile Prime Plus® VET, Cartucho automático de controles de química, de CO-oxímetro e de gás no sangue Stat Profile Prime Plus® VET com creatinina, Stat Profile Prime Plus® VET vérégáz, CO-oximéter, kémiai kontrollok automatikus patron kreatininél, (מחשנית עם קריאטינין) Stat Profile Prime Plus® VET Chemistry Controls Auto-Cartridge, Stat Profile Prime Plus® VET 血液ガス, CO オキシメーター, 生化学検査用コントロール自動カートリッジ (クレアチニン), 크레아티닌 사용 Stat Profile Prime Plus® VET 혈액 가스, CO-산소 농도계, 화학 조절제 자동 카트리지, Stat Profile Prime Plus® VET 血气、一氧化碳-血氧仪、化学对照溶液自动试剂盒 (含肌酐)

LOT 25211008

CONTROL 1 2 3 4 5

2027-01-14

Table with columns for Analyte, Unit, and Expected Ranges for Control 1, 2, 3, 4, and 5. Rows include pH, H+, PCO2, PO2, SO2, Hct, Na+, K+, Cl-, iCa, iMg, Glu, Lac, BUN, Urea, Creatinine, HbF, tHb, tHb, tHb, O2Hb, COHb, MetHb, HHb, tBil, and tBil.

EN

Product Description: Aqueous quality control material for monitoring the performance of pH, PCO2, PO2, SO2, hematocrit (Hct), fetal hemoglobin (HbF), total hemoglobin (tHb), oxyhemoglobin (O2Hb), carboxyhemoglobin (COHb), methemoglobin (MetHb), and deoxyhemoglobin (HHb) in Levels 1, 2 and 3 as well as Na+, K+, Cl-, iCa, iMg, Glucose, Lactate, BUN (urea), and Creatinine in Level 4 and 5. A charge kit consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation.

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Traceability of Standards: Total Hemoglobin (tHb) and Methemoglobin (MetHb) are traceable by using Cyanmethemoglobin method. Carboxyhemoglobin (COHb) and Oxyhemoglobin (O2Hb) are traceable using Spectrophotometry. Analytes are traced to NIST Standard Reference Materials.

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
1/NCCLS Document M29-T2. 2/How to Define and Determine Reference Intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13.

1/NCCLS Document M29-T2. 2/How to Define and Determine Reference Intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13.

1/NCCLS Document M29-T2. 2/How to Define and Determine Reference Intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13.

200 Prospect Street Waltham, MA 02454 U.S.A.

LPN 59633F 2024-10