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Evaluation of a Point-of-Care Glucose Meter for General Use in Complex Tertiary Care Facilities²⁰

Comparison of Four Hospital Based Glucose Meter Technologies for Accuracy, Precision and Interferences Encountered in Hospitalised Patients²¹

UCI neonatal

Performance of the Nova StatStrip Point-of-Care Glucose Meter in a Neonatal Intensive Care Unit²⁴

Clinical Performance of the New Glucometer in the Nursery and Neonatal Intensive Care Unit²⁵

UCI

Improved Blood Glucose Levels Achieved in ICU Patients Using Haematocrit Corrected Glucose Meter and Blood Gas Analyser Results²²

Validation of a Glucose Meter at an Intensive Care Unit²³

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Haematocrit Effects Lead to Inadequate Glycemic Control and Insulin Dosing in Adult Burn Patients²⁶

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Batería:

Tipo:..... 2 pilas AAA

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- https://trulst.co/blog/linux-statistics
- Diabetes UK. Diabetic ketoacidosis (DKA). 2016. https://www.diabetes.org.uk/guide-to-diabetes/complications/diabetic_ketoacidosis.
- Reaves M et al. ISPAD clinical practice consensus guidelines 2014 compendium. Assessment and monitoring of glycaemic control in children and adolescents with diabetes. *Pediatr Diabetes* 2014;15(Suppl 20):102-114.
- American Diabetes Association. Test of glycemia in diabetes. *Diabetes Care* 2004;17(Suppl 1):S91-93
- Vanelli M et al. Cost effectiveness of the direct measurement of 30-B-hydroxybutyrate in the management of diabetic ketoacidosis in children [Carta al editor]. *Diabetes Care* 2003;26:959.
- Noyes K et al. Hydroxybutyrate near-patient testing to evaluate a new end-point for intravenous insulin therapy in the treatment of diabetic ketoacidosis in children. *Pediatr Diabetes* 2007;150:156.
- Orjia IB et al. Pseudo hyperglycemia Secondary to High-Dose Intravenous Vitamin C Managed as Diabetic Ketoacidosis: An Endocrinological Catastrophe. *AAE Clin Case Reports* 2021.
- Krabbe J et al. Significant interference on specific point-of-care glucose measurements due to high dose of intravenous vitamin C therapy in critically patients. *Clin Chem Lab Med*. 59(5) 2021.
- Katzman B et al. Unintended Consequence of High-Dose Vitamin C Therapy for an Oncology Patient: Evaluation of Ascorbic Acid Interference with Three Hospital-Use Glucose Meters. *J Diabet Sci and Tech* 15(4) 2021.
- Lachance O et al. High-dose vitamin-C induced prolonged factitious hyperglycemia in a peritoneal dialysis patient: a case report. *J Med Case Reports* 2021.
- Dwan Ten Berge et al. Delayed Diagnosis of Severe Hypoglycemia in a Septic Patient With Chronic Renal Failure. *Cureus*. 14(8) 2022.
- Orjia IB et al. Pseudohyperglycemia Secondary to High-Dose Intravenous Vitamin C Managed as Diabetic Ketoacidosis: An Endocrinological Catastrophe. *AAE Clin Case Reports* 2021
- Pillay, T. S. et al. (2013) Continuing danger of glucose point-of-care test devices in the neonatal setting [Carta al editor]. *S Afr Med J*, 103(6), 356. doi:10.7196/SAMJ6918
- Kim, S. K. et al. (2013). Spurious elevation of glucose concentration during administration of high dose of ascorbic acid in a patient with type 2 diabetes on hemodialysis. *Tonsei Med J*, 54(5), 1289-1292. doi: 10.3349/jym.2013.54.12893.
- Mathew, V. et al. (2013). Erroneous glucose recordings while using mutant variant of quinoprotein glucose dehydrogenase glucometer in a child with galactosemia. *Indian J Endocr Metab*, 17(Suppl. 1), 289-291. doi:10.4103/2230-8210.119616
- Harper C. FDA perspective - Public Health Notification: Potentially fatal errors with GDH-POD glucose monitoring technology. Presented at the FDA/CDRH public meeting, Gaithersburg, MD: 2010.
- Darragh T. Report: Lehigh Valley Hospital gave fatal dose. *The Morning Call*, 2011 Jun 15.
- Frias J et al. Review of adverse events associated with false glucose readings measured by GDH-POD-based glucose test strips in the presence of interfering sugars. *Diabetes Care* 2010;33:728-729.
- Food and Drug Administration. FDA clears glucose monitoring system for use in hospital critical care units [Comunicado de prensa]. 2014 Sep 24.
- Chan P et al. Evaluation of a point-of-care glucose meter for general use in complex tertiary care facilities. *Clin Biochem* 2009;42:1104-1112.
- Bewley B et al. (Junio de 2007). Comparison of four hospital based glucose meter technologies for accuracy, precision and interferences encountered in hospitalized patients. Sesión con láminas presentada en la reunión del European Congress of Clinical Chemistry and Laboratory Medicine, Amsterdam, Países Bajos.
- Roman A et al. (Febrero de 2011). Improved blood glucose levels achieved in ICU patients using hematocrit corrected glucose meter and blood gas analyzer results. Sesión con láminas presentada en la reunión de Advanced Technologies & Treatments for Diabetes, Londres, Reino Unido.
- López M et al. Validation of a glucose meter at an intensive care unit. *Endocrinol Nutr* 2012;59:28-34.
- Tendil K et al. (Octubre de 2012). Performance of the Nova StatStrip point of care glucose meter in a neonatal intensive care unit. Sesión con láminas presentada en la reunión del European Joint Congress, Dubrovnik, Croacia.
- Nunamaram P et al. Clinical performance of the new glucometer in the nursery and neonatal intensive care unit. *Pediatr Int* 2011; 53:218-223.
- Godwin Z et al. (Julio de 2012). Hematocrit effects leads to inadequate glycaemic control and insulin dosing in adult burn patients. Sesión con láminas presentada en la reunión de la Asociación Americana de Química Clínica, Los Angeles, California.
- Tran N et al. Clinical impact of sample interference on intensive insulin therapy in severely burned patients: A pilot study. *J Burn Care Res* 2014;35:72-79.
- Boyd J et al. Quality specifications for glucose meters: Assessment by simulation modeling of errors in insulin dose. *Clin Chem* 2001;47:209-214.
- Hoff S et al. Comparison of point-of-care testing glucose results from intensive care patients measured with network-ready devices. *Diabetes Technol Ther* 2011;13:1047-1056.
- Mann E et al. Error rates resulting from anemia can be corrected in multiple commonly used point-of-care glucometers. *J Trauma* 2006;64:15-20.
- Corwin H et al. The CRIT Study: Anemia and blood transfusion in the critically ill—Current clinical practice in the United States. *Crit Care Med* 2004;32:39-52.
- Karon B et al. Evaluation of the impact of hematocrit and other interference on the accuracy of hospital-based glucose meters. *Diabetes Technol Ther* 2008;10:111-120.
- Mohn B. Evaluation of the Nova Biomedical StatStrip glucose meter. *NZ J Med Lab Sci* 2010;64:18-21.
- Raizman J et al. Impact of improved glucose monitoring in the neonatal intensive care unit: An evaluation of the clinical performance of the point of care Nova StatStrip glucose meter. Límina presentada en la Reunión anual de AAC, Chicago, IL: 2014.
- Raizman J et al. Clinical impact of improved point-of-care glucose monitoring in neonatal intensive care using Nova StatStrip: Evidence for improved accuracy, better sensitivity, and reduced test utilization. *Clin Biochem* 2016;49:879-884.
- Karon B et al. Impact of glucose meter error on glycemic variability and time in target range during glycaemic control after cardiovascular surgery. *J Diabetes Sci Technol* 2016;10:336-342.
- Karon B et al. Retrospective evaluation of the accuracy of Roche AccuCheck Inform and Nova StatStrip glucose meters when used on critically ill patients. *Diabetes Technol Ther* 2014;16:1-5.
- DuBois J et al. Bedside glucose monitoring—is it safe? A new, regulatory-compliant risk assessment evaluation protocol in critically ill patient care settings. *Crit Care Med* 2017;45:567-574.
- Food and Drug Administration. Blood glucose monitoring test systems for prescription point-of-care use. Guidance for industry and Food and Drug Administration staff. Silver Spring, MD: FDA, 2016.
- DuBois J et al. Multi-site evaluation of a hospital POC blood glucose monitoring system in critically ill patients on intensive insulin therapy. Límina presentada en la reunión internacional Hospital Diabetes Meeting, Atlanta, GA.
- Baum J et al. Improving the quality of self-monitoring blood glucose measurement: A study in reducing calibration errors. *Diabetes Technol Ther* 2006;8:237-247.

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nova
 biomedical
 novabiomedical.com

Nova Biomedical Headquarters: 200 Prospect St., Waltham, MA 02454 U.S.A., +1-781-894-0800 800-458-5813, FAX: +1-781- 894-5915, Int'l FAX: +1-781-899-0417, e-mail: info@novabio.com
Nova Biomedical Australia ANZ Pty. Ltd.: 5/372 Eastern Valley Way, Chatswood, NSW, 2067, Australia, TEL: +61(0) 2 9417 0193, e-mail: AU-info@novabio.com
Nova Biomedical Benelux B.V.: Korenolen 22, 5281 PB, Boxtel, The Netherlands, TEL: +31(0)733032701, e-mail: benelux-info@novabio.com
Nova Biomedical Brasil: Rua Massena, 107, Jardim Canadá, Nova Lima - MG, CEP: 34007-746 Brasil, TEL: +55-31-3360-2500, e-mail: BR-info@novabio.com
Nova Biomedical Canada, Ltd: 17 - 2900 Argenta Road, Mississauga, Ontario L5N 7X9 Canada, TEL: +1-905-567-7700 800-263-5999, FAX: +1-905-567-5496, e-mail: CA-info@novabio.com
Nova Biomedical France: Parc Technopolis - Bât. Sigma 3 Avenue du Canada 91940 Les Ulis Courtabouef, France, TEL: +33-1-64 86 11 74, FAX: +33-1-64 46 24 03, e-mail: FR-info@novabio.com
Nova Biomedical GmbH, Deutschland: Hessenring 13 A, Geb. G, 64546 Mörfelden-Walldorf, Germany, TEL: +49-6105 4505-0, FAX: +49-6105 4505-37, e-mail: DE-info@novabio.com
Nova Biomedical Iberia, S.L.: c/Vic 17, Planta 3A 08173 Sant Cugat del Vallés, Barcelona, Spain, TEL: +34 935531173, e-mail: ES-info@novabio.com or PT-info@novabio.com
Nova Biomedical Italia S.r.l.: via Como, 19 - 20045 Lainate (MI), Italy, Tel: +39 02 87070041, Fax: +39 02 87071482, e-mail: IT-info@novabio.com
Nova Biomedical K.K., Japan: Harumi Island Triton Square Office Tower X 7F, 1-8-10 Harumi, Chuo-ku, Tokyo 104-6007, Japan, TEL: 03-5144-4144, FAX: 03-5144-4177, e-mail: jp-info@novabio.com
Nova Biomedical New Zealand: Regus Constellation Drive - Cantida Building 4 Level 3/61 Constellation Drive, Rosedale, Auckland, 0630, New Zealand, TEL: +64 800 555 268, e-mail: nz-sales@novabio.com
Nova Biomedical Schweiz GmbH: Herostrasse 7, 8048 Zürich, Switzerland, TEL: +41-41-521-6655, FAX: +41-41-521-6656, e-mail: CH-info@novabio.com
Nova Biomedical U.K.: Innovation House, Aston Lane South, Runcorn, Cheshire WA7 3FY United Kingdom, TEL: +44-1928 704040, FAX: +44-1928 796792, e-mail: UK-info@novabio.com

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