# Stat Profile® pHox Ultra/Critical Care Xpress

## Chemistry Controls & Performance Check Solution Auto-Cartridge

Chemiekontrollen und Leistungstest-Lösung Auto-QK-Kassette, Αυτόματη Κασέτα Χημικών Ελεγκτών & Διαλύματος Ελέγχου Απόδοσης, Cartucho automático para Controles para química y Solución de verificación de comportamiento, Cartouche CQ de contrôles chimiques et de solution de contrôle de linéarité, Cartuccia automatica controlli chimici e soluzione controllo prestazioni, Controlos de química & Cartucho automático de Solução de Verificação de Desempenho, Autokassett med kemisk kontroller och prestandakontrollösnig

LOT 25302061

2027-04-21

CONTROL 4 5

PERFORMANCE CHECK

SOLN

### Expected Ranges, Erwartete Bereiche, Αναμενόμενα Εύρη, Rangos esperados, Fourchettes attendues, Intervalli previsti, Gamas previstas, Förväntade områden

		CONTROL 4		CONTROL 5			
		min	$\bar{x}$	max	min	$\bar{x}$	max
Na <sup>+</sup>	mmol/L	139.2	143.2	147.2	112.8	116.8	120.8
K <sup>+</sup>	mmol/L	3.80	4.05	4.30	6.28	6.58	6.88
Cl <sup>-</sup>	mmol/L	122.2	126.7	131.2	91.5	96.0	100.5
Ca <sup>++</sup>	mmol/L	1.03	1.11	1.19	1.52	1.62	1.76
Mg <sup>++</sup>	mmol/L	0.54	0.63	0.72	0.98	1.13	1.28
Glu	mg/dL	73	80	87	255	280	305
Glu	mmol/L	4.05	4.44	4.83	14.15	15.54	16.93
Lac	mmol/L	1.7	2.0	2.3	6.3	7.0	7.7
BUN	mg/dL	14	19	24	49	56	63
BUN	mmol/L	5.0	6.8	8.6	17.5	20.0	22.5
Urea	mg/dl	30.0	40.8	51.5	105.1	120.1	135.1
Urea	mmol/L	5.0	6.8	8.6	17.5	20.0	22.5

PERFORMANCE CHECK	SOLN
GLUCOSE	PASS - No Action
	FAIL - Replace Membrane
	FEHLERFREI - Keine Maßnahme
	FEHLERHAFT - Membran austauschen
	ΕΠΙΤΥΧΙΑ - Καμία Ενέργεια
	ΑΠΟΤΥΧΙΑ - Αντικαταστήστε τη μεμβράνη
	PASÓ - Ninguna acción
	NO PASÓ - Sustituir membrana
	SUCCES - Aucune action
	ECHEC - Remplacer la membrane
RIUSCITO - Nessuna azione	
NON RIUSCITO - Sostituisci membrana	
PASSOU - Nenhuma acción	
FALHOU - Substituir membrana	
OK - Ingen åtgärd	
FEL - Byt ut membran	

#### Product Description

Carton contains three bags. Two bags contain assayed quality control material formulated for monitoring the measurement of Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, Ca<sup>++</sup>, Mg<sup>++</sup>, Glucose, Lactate, and BUN (Controls 4 and 5). One bag contains an aqueous quality control material to verify glucose membrane performance (Performance Check Solution, Control 6).

#### Intended Use

For in vitro diagnostic use to monitor the performance of the Stat Profile pHox Ultra/Critical Care Xpress Analyzer and to insure optimum performance of glucose membranes. Performance check solution should be analyzed before controls are run and after changing glucose membrane.

#### Methodology

Refer to Stat Profile pHox Ultra/Critical Care Xpress Analyzer Instruction for Use Manual.

#### Composition

Controls (Controls 4 and 5) contain known concentrations of Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, Ca<sup>++</sup>, Mg<sup>++</sup>, Glucose, Lactate, and BUN at two levels. Performance Check Solution (Control 6) is a buffered salt solution containing Glucose (100mg/dL). Each bag contains a minimum volume of 90 mL. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials.

#### Warning and Precautions

For in vitro diagnostic use. Do not freeze. Refer to the Stat Profile pHox Ultra/Critical Care Xpress Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory controls.

#### Storage

Store at 2-8°C. Do not freeze. Each cartridge has lot number and expiration date printed on the label.

#### Directions for Use

Refer to Analyzer Instructions for Use Manual for complete directions. Performance Check solution should be analyzed before controls are run and after changing glucose membranes. Verify that the lot number appearing on the expected ranges table is the same as the lot number printed on the cartridge.

#### Limitations

For use on Nova Biomedical Analyzers Only. Performance characteristics have not been established for use on other manufacturers' analyzers. Once installed, each Control Auto-Cartridge may be used for a maximum of 35 days from the date of initial activation and installation. The Analyzer will indicate that the cartridge is invalid at that time.

This product is warranted to perform as described in this product insert, and Nova Biomedical disclaims any implied warranty of merchantability or fitness for any other purpose, and in no event shall Nova Biomedical be liable for any consequential damages arising out of the aforesaid expressed warranty.

#### Traceability of Standards

Analytes are traced to NIST Standard reference materials.

#### Reference Intervals

Concentrations are formulated at normal and abnormal levels expected in patient blood. The expected clinical ranges of patient blood is referenced in Tietz, N.W. ed. 1986 Textbook of Clinical Chemistry, W.B. Saunders Co. Benutzer wish to determine Expected Ranges in their own laboratory.<sup>2</sup>

#### Expected Ranges

The expected range for each analyte was determined at Nova Biomedical using replicate determinations on Nova analyzers. The expected range indicates the maximum deviations from the mean value that may be expected under differing laboratory conditions for analyzers operating within specifications. Refer to Expected Ranges table for Control 4 and 5 values. Performance Check Solution results for glucose will be Pass or Fail.

<sup>1</sup> NCCLS Document M29-T2.

<sup>2</sup> How to define and determine reference intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13.

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