REF 48559, 41549

IVD 15°C 30°C CE

Stat Profile[®] pHox Ultra/Critical Care Xpress CO-Oximeter Calibrator Cartridge with Bilirubin and Deproteinizing Solution

Kalibrator-Patrone mit Bilirubin und Deproteinisierungslösung, Πακέτο βαθμονομητή με χοληρυθρίνη και διάλυμα αποπρωτεϊνισμού, Cartucho de calibrador con bilirubina y solución desproteinizante, Cartouche de calibrateur avec bilirubine et solution de déprotéinisation, Cartuccia calibratore con bilirubina e soluzione deproteinizzante, Cartucho de calibrador com solução de desproteinização e bilirubina, ビリルビンおよび除蛋白液付きキャリブレータ・カートリッジ, Calibrator Cartridge와 빌리루빈 및 단백질 제거액, 含有胆红素和去蛋白溶液的定标品盒

LOT 25321003 2026-11-17

tHb CALIBRATOR ≥70mL		DIL ≥180mL	SOLN F ≥620mL	Deproteinizing Solution ≥90mL	WASTE CONTAINER W						
Units <table border="1"> <tr> <td>g/dL</td> <td>16.2</td> </tr> <tr> <td>g/L</td> <td>162</td> </tr> <tr> <td>mmol/L</td> <td>10.0</td> </tr> </table>		g/dL	16.2	g/L		162	mmol/L	10.0	<p>Deionized Water, Surfactant, Buffers, Preservatives, Deionat, Tensid, Puffer, Konservierungsmittel Αποιονισμένο νερό, επιφανειοδραστικό μέσο, ρυθμιστικά διαλύματα, συντηρητικά, Agua desionizada, surfactante, tampones, conservantes, Eau désionisée, agent de surface, tampons, agents de conservation, Acqua deionizzata, surfattante, tamponi, conservanti, Água desionizada, Surfactante, Tampões, Conservantes, 脱イオン水、界面活性剤、緩衝剤、保存剤 탈이온수, 계면 활성제, 버퍼, 방부제 去离子水、表面活性剂、缓冲剂、防腐剂</p> <p>Deionized Water, Surfactant, Preservatives, Deionat, Tensid, Konservierungsmittel, Αποιονισμένο νερό, επιφανειοδραστικό μέσο, συντηρητικά, Acqua desionizzata, surfactante, conservantes, Eau désionisée, agent de surface, agents de conservation, Acqua deionizzata, surfattante, conservanti, Água desionizada, Surfactante, Conservantes, 脱イオン水、界面活性剤、保存剤 탈이온수, 계면 활성제, 방부제 去离子水、表面活性剂、防腐剂</p> <p>Deproteinizing Solution Deproteinisierungslösung Διάλυμα αποπρωτεϊνισμού Solución desproteinizante Solution de déprotéinisation Soluzione deproteinizzante Solução de desproteinização 除蛋白液 단백질 제거액 去蛋白溶液 NaOCl 0.4 mol/L</p>		
g/dL	16.2										
g/L	162										
mmol/L	10.0										

Description
 For use in calibrating the CO-Oximeter on the Stat Profile pHox Ultra/Critical Care Xpress (CCX) Analyzer. The above table lists the concentrations of the constituents for the calibrator, flush, deproteinizing and diluent solutions contained in the Calibrator Cartridge. Contains the following aqueous solutions within individual flexible bags housed in a cardboard box: flush solution, diluent solution, deproteinizing solution, calibrator and a flexible waste bag. The exposed bag fittings are arranged in a line along the rear of the cartridge box. The septa are pierced during the insertion of the cartridge into the analyzer.

Intended Use
 For the quantitative determination of total hemoglobin, oxyhemoglobin, carboxyhemoglobin, methemoglobin, deoxyhemoglobin and total bilirubin in human blood using the Stat Profile pHox Ultra/Critical Care Xpress CO-Oximeter Analyzer.

Methodology
 Refer to the Stat Profile pHox Ultra/Critical Care Xpress Analyzer Instructions for Use Manual for Test Methodologies and Principles.

Composition
 Solution volumes are listed in the table above. Contains no constituents of human origin, however good laboratory practice should be followed during handling of these materials. (REF: NCCLS DOCUMENT M29-T2).

Warnings and Cautions
 Intended for in vitro diagnostic use. Refer to the Stat Profile pHox Ultra/Critical Care Xpress Analyzer Instructions for Use Manual for complete directions. Follow standard practices for handling laboratory reagents. Discard Calibrator Cartridges in accordance with local guidelines. Waste bag contains Preservative (eye, skin and respiratory irritant).

Storage
 Store at 15°C - 30°C. Expiration date is printed on the front of the cartridge.

Directions for Use
 Refer to the Stat Profile pHox Ultra/Critical Care Xpress Analyzer Instructions for Use Manual for complete directions for installation of the cartridge. Cartridge is ready for use.

Limitations
 For use on Nova Biomedical Analyzers only.
 If not stored at room temperature and out of direct sunlight, the solution may deteriorate. Indications of deterioration may be high or low out of range control values or a noticeable change in calibration or both. Cloudiness or particulate matter may be a sign of possible deterioration.
 This product is warranted to perform as described in this product insert, and Nova Biomedical disclaims any implied warranty or merchantability or fitness for any other purpose, and in no event shall Nova Biomedical be liable for any consequential damages arising out of the aforesaid expressed warranty.

Traceability of Standards
 Traceable to precise hemoglobin measurement using the reference cyanmethemoglobin method.

Reference
 ICSH, Recommendations for Haemoglobinometry in Human Blood. British Journal of Haematology, 13 (suppl.) pp 71-75, 1967.

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