Stat Profile Prime Plus® Chemistry Controls 4 and 5

Controles de química 4 y 5 Stat Profile Prime Plus®, Contrôles 4 et 5 chimie Stat Profile Prime Plus®, Stat Profile Prime Plus® Blutchemiekontrollen 4 und 5, Υλικό ελέγχου Χημείας 4 και 5 Stat Profile Prime Plus®, Controlli chimici 4 e 5 Stat Profile Prime Plus®, Controles de Química Stat Profile Prime Plus® 4 e 5, Stat Profile Prime Plus® Kimya Kontrolleri 4 ve 5, Stat Profile Prime Plus® 4-es és 5-ös kémiai kontroll, 5 - 1 4 בקרת כימיה Stat Profile Prime Plus®, Stat Profile Prime Plus® 生化学検査用コントロール 4 および 5, Stat Profile Prime Plus® 화학 조절제 4 및 5, Stat Profile Prime Plus® 化学对照溶液 4 和 5

CONTROL 4 5

LOT 25227026 2026-08-13

EN

Product Description
Aqueous quality control material for monitoring the performance of Na⁺, K⁺, Cl⁻, iCa, iMg, Glucose, Lactate, BUN (urea), and Creatinine for use with Stat Profile Prime Plus Analyzers ONLY.
Formulated at 2 levels

Intended Use
Intended for *in vitro* diagnostic use by healthcare professionals use for monitoring the performance of the Stat Profile Prime Plus Analyzers.

Methodology
Refer to Stat Profile Prime Plus Analyzer Instruction for Use Manual for Methodology and Principles.

Composition
Controls are buffered solutions containing known concentrations of Na⁺, K⁺, Cl⁻, iCa, iMg, Glucose, Lactate, BUN (Urea), Creatinine, and preservatives. Each ampule contains a minimum volume of 1.7 mL. Controls contain no constituents of human origin, however good laboratory practices should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-T21¹).

Warning and Precautions
Controls must be stored at 2-8°C (37-46°F). Once opened, analyze within one minute and discard the unused portion in accordance with local guidelines. Refer to Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete information. Intended for *in vitro* diagnostic use. Follow standard practices for handling laboratory reagents.

Storage
Store at 2-8°C (37-46°F). DO NOT FREEZE. Each ampule has a Lot Number and Expiration Date printed on the label.

Directions for Use
Store controls at 2-8°C (37-46°F) until ready for use. Ensure control is at room temperature prior to use. Shake ampule well before opening, snap open ampule (protecting fingers with gauze or glove). Once opened, analyze control within one minute. Discard the unused portion in accordance with local guidelines. Verify that the Lot Number on the Expected Ranges Table corresponds to the Lot Number on the ampule. Refer to Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete instructions.

Limitations
The Expected Range values are specific for instruments and controls manufactured by Nova Biomedical.

Traceability of Standards
Analytes are traced to NIST Standard reference materials.

Reference Intervals
Concentrations are formulated at normal and abnormal expected values in patient blood. The expected clinical range of these values in patient blood is referenced in Tietz, NW ed. 1986 Textbook of Clinical Chemistry, W.B. Saunders Co. Users may wish to determine Mean Values and Expected Ranges in their own laboratory.²

Expected Ranges
The expected range for each parameter was determined at Nova Biomedical using replicate determinations on Nova analyzers. The expected range indicates the maximum deviations from the Mean Value that may be expected under differing laboratory conditions for instruments operating within specifications. Refer to Expected Ranges Table.

¹NCCLS Document M29-T2.
²How to Define and Determine Reference Intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13.

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