REF 49187

IVD 15°C 30°C

Stat Profile® pHox Ultra/Critical Care Xpress CO-Oximeter Controls

CO-Oximeter-Kontrollen, Διαλύματα Ελέγχου Οξυμέτρου-CO, Controles de CO-Oxímetro, Contrôles par oxymètre CO, Controlli qualità co-ossimetro, Controles de co-oxímetro, CO-oximeter-kontroller, CO-酸素濃度計コントロール, CO-산소 농도계 조절, CO-血氧計质控品

CONTROL 7 8 9

LOT 22188032

2024-06-24

Table with columns for LOT, CONTROL (7, 8, 9), and Expected Ranges for parameters like tHb, O2Hb, MetHb, HHb, tBil. Includes sub-headers in multiple languages: Expected Ranges, Erwartete Bereiche, Αναμενόμενα εύρη, Rangos esperados, Fourchettes attendues, Intervalli previsti, Gamas previstas, Förväntade områden, 予測範囲, 예상 범위, 预期范围.

Description: An aqueous quality control material for monitoring the measurement of total hemoglobin (tHb), oxyhemoglobin (O2Hb), carboxyhemoglobin (COHb), methemoglobin (MetHb), deoxyhemoglobin (HHb) and total bilirubin (tBil).

Intended Use: For in vitro diagnostic use by health care professionals for monitoring the performance of Stat Profile pHox Ultra/CCX Analyzers.

Methodology: Refer to Stat Profile pHox Ultra/CCX Analyzer's Instructions for Use Manual for use, Methodology and Principles.

Composition: Contains dyes, viscosity enhancers and a preservative in an aqueous salt solution. Each ampule contains 1.7 mL. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials.

Warnings and Cautions: Intended for in vitro diagnostic use. Refer to the Stat Profile pHox Ultra/CCX Analyzer's Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory reagents.

Storage: Store at 15-30°C. Exposure to direct sunlight should be avoided. Expiration date is printed on each ampule. Do not store unused opened ampules.

Directions for use: Refer to the analyzer's Instructions for Use Manual for complete directions. Each ampule is intended for single use only. However, multiple analyses can be processed from one ampule.

Limitations: This product is warranted to perform as described in this product insert, and Nova Biomedical disclaims any implied warranty or merchantability or fitness for any other purpose.

Traceability of Standards: Total Hemoglobin (tHb) and Methemoglobin (MetHb) are traceable by using Cyanmethemoglobin method. Carboxyhemoglobin (COHb) and Oxyhemoglobin (O2Hb) are traceable using Spectrophotometry.

Reference Intervals: Concentrations are formulated to represent three levels.

Users may wish to determine MEAN VALUES and EXPECTED RANGES in their own laboratory*.

Expected Ranges: The EXPECTED RANGE for each parameter was determined at Nova Biomedical by using multiple runs of each level of control on properly calibrated instruments.

*How to Define and Determine Reference Intervals in the clinical laboratory; approved guideline-second edition, NCCLS C28-A2, Volume 20, Number 13

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