REF 52427,52861,53365,53105,53469,53470,52862,52863,53364,53466,53467,53468

Rx only IVD 15°C 30°C CE

# Stat Profile® Prime Calibrator Cartridge CCS/CCS Comp

Kalibratorkassette, Κασέτα Βαθμονομητή, Cartucho del calibrador, Cartouche étalon, Cartuccia calibratore, Calibrator Cartridge, Kalibratorkassett, Kalibratorkazetta, キャリブレータ カートリッジ, 측정기 카트리지, 校准液盒

LOT 23060084

2024-06-08

CAL A ≥ 500 mL	CAL B ≥ 250 mL	CAL F ≥ 720 mL	SOLN R ≥ 300 mL
pH 7.360	pH 6.840	PCO <sub>2</sub> 51.5 mmHg	KCl 2 mol/L
PCO <sub>2</sub> 24.3 mmHg	PO <sub>2</sub> 167.0 mmHg	PO <sub>2</sub> 93.9 mmHg	
Na <sup>+</sup> 133.0 mmol/L	Na <sup>+</sup> 68.5 mmol/L		
K <sup>+</sup> 4.0 mmol/L	K <sup>+</sup> 10.0 mmol/L		
Cl <sup>-</sup> 106.0 mmol/L	Cl <sup>-</sup> 46.0 mmol/L		
iCa 1.10 mmol/L	iCa 2.20 mmol/L		
Glu 80 mg/dL	Glu 200 mg/dL		
Lac 2.0 mmol/L	Lac 10.0 mmol/L		

(EN)

**Product Description**

The above table lists the concentrations of constituents for the calibrators and solutions contained in the Calibrator Cartridge. Contains the following aqueous solutions within individual flexible bags held in a cardboard box: 3 Calibrators (Calibrator A, B, and F), reference materials and a flexible waste bag. Each bag includes a filament with a septa. The exposed bag filaments are arranged in a line along the rear of the cartridge box. These septa are pierced during the insertion of the cartridge into the analyzer.

**Intended Use**

The Stat Profile Prime Calibrator Cartridge CCS is intended for the calibration of pH, PCO<sub>2</sub>, PO<sub>2</sub>, Hct, Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, iCa, Glucose (Glu), and Lactate (Lac) using the Stat Profile Prime

**Methodology**

Refer to Prime Analyzer Instructions For Use Manual for Test Methodologies and Principles.

**Composition**

Calibrator volumes are listed in the table above. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-T2).

**Warnings and Cautions:**

Intended for in vitro diagnostic use.  
Refer to the Prime Analyzer Instructions for Use Manual for complete directions for use, installation instructions, recommendations for calibration intervals, Troubleshooting information, Methodology and Principles of testing procedures. Follow standard practices required for handling laboratory reagents. Discard Calibrator Cartridges and Waste products in accordance with local guide lines. Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

**Storage**

Store at 15-30°C. Expiration date is printed on the front of the Calibrator Cartridge.

Directions For Use Refer to the Prime Analyzer Instructions for Use for complete directions for installation of the cartridge, and recommendations for Calibration Intervals. Calibrator Cartridges must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use.

**Traceability of Standards**

Analyses traced to NIST Standard Reference Materials.

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