REF 53106

Stat Profile® Prime VET CCS Comp Calibrator Cartridge

Kalibratorkassette, Κάσета Βαθμονομητή, Cartucho del calibrador, Cartouche étalon, Cartuccia calibratore, Calibrator Cartridge, Kalibratorkassett, Kalibratorkazetta, 교정기 카트리지, キャリブレータカートリッジ, 校准液盒

LOT 23060082

2024-06-08

CAL A ≥ 500 mL		CAL B ≥ 250 mL		CAL F ≥ 720 mL		SOLN R ≥ 300 mL	
pH	7.360	pH	6.840	PCO ₂	51.5 mmHg	KCl	2 mol/L
PCO ₂	24.3 mmHg	PO ₂	167.0 mmHg	PO ₂	93.9 mmHg	 WASTE CONTAINER W	
Na ⁺	133.0 mmol/L	Na ⁺	68.5 mmol/L				
K ⁺	4.0 mmol/L	K ⁺	10.0 mmol/L				
Cl ⁻	106.0 mmol/L	Cl ⁻	46.0 mmol/L				
iCa	1.10 mmol/L	iCa	2.20 mmol/L				
Glu	80 mg/dL	Glu	200 mg/dL				
Lac	2.0 mmol/L	Lac	10.0 mmol/L				

EN
Product Description
 The above table lists the concentrations of constituents for the calibrators and solutions contained in the Calibrator Cartridge. Contains the following aqueous solutions within individual flexible bags housed in a cardboard box: 3 Calibrators (Calibrator A, B, and F), reference solution, and a flexible waste bag. Each bag includes a filament with a septa. The exposed bag filaments are arranged in a line along the rear of the cartridge box. These septa are pierced during the insertion of the cartridge into the analyzer.

Intended Use
 For the quantitative determination of Na, K, Cl, iCa, Glucose, Lactate, pH, PCO₂, PO₂, and Hct in blood using the Stat Profile Prime VET CCS Comp Analyzer.

Methodology
 Refer to Prime Analyzer Instructions For Use Manual for Test Methodologies and Principles.

Composition
 Calibrator volumes are listed in the table above. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF. NCCLS DOCUMENT M29-T2).

Warnings and Cautions:
 Refer to the Prime Analyzer Instructions for Use Manual for complete directions for use, installation instructions, recommendations for calibration intervals, troubleshooting information, Methodology and Principles of the testing procedures. Follow standard practices required for handling laboratory reagents. Discard Calibrator Cartridges and Waste products in accordance with local guide lines. Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

Storage
 Store at 15-30°C. Expiration date is printed on the front of the Calibrator Cartridge.

Directions For Use
 Refer to the Prime Analyzer Instructions for Use for complete directions for installation of the cartridge, and recommendations for Calibration Intervals. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use.

Traceability of Standards
 Analytes traced to NIST Standard Reference Materials.

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