REF 57825, 57826, 57827, 57828, 57829

RX Only IVD 15°C 30°C CE

Stat Profile Prime Plus® Calibrator Cartridge

Stat Profile Prime Plus® Kalibrierkassette, Φυσιγγιο μέσου βαθμονόμησης με κρεατινίνη Stat Profile Prime Plus®, Cartouche de calibrage Stat Profile Prime Plus®, Cartouche de calibrage Stat Profile Prime Plus®, Cartuccia calibratore Stat Profile Prime Plus®, Cartucho de calibrador Stat Profile Prime Plus®, Stat Profile Prime Plus®' kalibrátor patron, כריכת בדיקת Stat Profile Prime Plus®, Stat Profile Prime Plus® キャリブレーターカートリッジ, Stat Profile Prime Plus® 보정기 카트리지, Stat Profile Prime Plus® 校准试剂盒

LOT 23089034

2024-08-21

CAL A	$\geq 495 \text{ mL}$	CAL B	$\geq 250 \text{ mL}$	CAL C	$\geq 245 \text{ mL}$	CAL D	$\geq 250 \text{ mL}$	CAL F	$\geq 850 \text{ mL}$	SOLN R	$\geq 500 \text{ mL}$
HCO ₃	11 mmol/L	pH	6.834 @ 37°C	HCO ₃	22 mmol/L	pH	7.354 @ 37°C	PO ₂	106.1 mmHg	KCl	2 mol/L
Na ⁺	135.5 mmol/L	H ⁺	146.55 nmol/L	Cl ⁻	69.0 mmol/L	H ⁺	44.26 nmol/L	PO ₂	14.11 kPa		
K ⁺	4.0 mmol/L	Na ⁺	77.0 mmol/L	Glu	200 mg/dL						
Cl ⁻	108.0 mmol/L	K ⁺	10.0 mmol/L	Glu	11.10 mmol/L						
iCa	1.04 mmol/L	iCa	1.92 mmol/L	Lac	10.0 mmol/L						
iCa	4.16 mg/dL	iCa	7.68 mg/dL	Lac	89.0 mg/dL						
iMg	0.50 mmol/L	iMg	1.50 mmol/L	BUN	50 mg/dL						
iMg	1.22 mg/dL	iMg	3.63 mg/dL	BUN	17.9 mmol/L						
Glu	80 mg/dL	tHb	15.2 g/dL								
Glu	4.44 mmol/L	tHb	152 g/L								
Lac	2.0 mmol/L	tHb	9.4 mmol/L								
Lac	17.8 mg/dL										
BUN	10 mg/dL										
BUN	3.6 mmol/L										

EN**Description**

The Cartridge contains the following aqueous solutions within individual bags: 5 Calibrators (Cartridges A, B, C, D and F) and one Preservative Bag. Each bag includes a septum. Each bag is supplied with a septum. The exposed bag filaments are arranged in a slot along the rear of the cartridge box. The septa are pierced during the insertion of the cartridge into the analyzer.

Intended Use

For the quantitative determination of Na⁺, K⁺, Cl⁻, iCa, iMg, BUN (Urea), Glucose, Lactate, pH, PCO₂, PO₂, SO₄²⁻, Hct, Hb, HbF, O₂Hb, COHb, MetHb, HHb, and iBil⁺ in human blood using the Stat Profile Prime Plus Analyzer.

*Not available in the USA or for Point-of-Care/Near-Patient Testing use.

Methodology

Refer to Stat Profile Prime Plus Analyzer Instructions For Use Manual for Methodologies and Principles.

Composition

Calibration and solution volumes are listed in the table. Contains no constituents of human origin, however good laboratory practice should be followed during handling of these materials. (REF NCCLS DOCUMENT M29-T2).

Warnings and Cautions**Intended for in vitro diagnostic use**

Refer to the Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory reagents. Discard Cartridge in accordance with local guidelines.

Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

Storage

Store at 15 - 30 °C. Expiration date is printed on the Calibrator Cartridge.

Directions For Use

Refer to the Stat Profile Prime Plus Analyzer Instructions for Use for complete directions.

Calibrator Cartridge must be inverted gently several times to thoroughly mix contents.

Traceability of Standards

Analytes traced to NIST Standard Reference Materials wherever possible.

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EC REP Nova Biomedical GmbH Hessenring 13A, Geb. G 64546 Mörfelden-Walldorf Germany
Waltham, MA 02454 U.S.A.