Stat Profile Prime Plus® VET Calibrator Cartridge With Creatinine

Stat Profile Prime Plus® VET Kalibrierkassette mit Kreatinin, Φυσίγγιο μέσου βαθμονόμησης VET με κρεατινίνη Stat Profile Prime Plus®, Cartucho de calibrador con creatinina Stat Profile Prime Plus® VET, Cartouche de calibrateurs avec créatinine Stat Profile Prime Plus® VET, Cartuccia calibratore con creatinina Stat Profile Prime Plus® VET, Cartucho de calibrador Stat Profile Prime Plus® VET com creatinina, Stat Profile Prime Plus® VET kalibrátor patron kreatininan, Stat Profile Prime Plus® VET מבחן ניטרול ציון קייל עם קריינטינין, Stat Profile Prime Plus@VET キャリブレーター(クリアチニン), 크레이티닌 사용 Stat Profile Prime Plus® VET 보정기 카트리지, Stat Profile Prime Plus® 含肌酐的校准试剂盒

LOT

23018041

2024-06-15

CAL A ≥ 495 mL	CAL B ≥ 250 mL	CAL C ≥ 245 mL	CAL D ≥ 250 mL	CAL F ≥ 850 mL	SOLN R ≥ 500 mL
HCO ₃ 11 mmol/L	pH 6.834 @ 37°C	HCO ₃ 22 mmol/L	pH 7.354 @ 37°C	Po ₂ 107.8 mmHg	KCl 2 mol/L
Na ⁺ 135.5 mmol/L	H ⁺ 146.55 nmol/L	Cl ⁻ 69.0 mmol/L	H ⁺ 44.26 nmol/L	Po ₂ 14.34 kPa	WASTE CONTAINER W
K ⁺ 4.0 mmol/L	Na ⁺ 77.0 mmol/L	Glu 200 mg/dL			
Cl ⁻ 108.0 mmol/L	K ⁺ 10.0 mmol/L	Glu 11.10 mmol/L			
iCa 1.04 mmol/L	iCa 1.92 mmol/L	Lac 10.0 mmol/L			
iCa 4.16 mg/dL	iCa 7.68 mg/dL	Lac 89.0 mg/dL			
iMg 0.50 mmol/L	iMg 1.50 mmol/L	BUN 50 mg/dL			
iMg 1.22 mg/dL	iMg 3.63 mg/dL	BUN 17.9 mmol/L			
Glu 80 mg/dL	tHb 15.2 g/dL	Creat 5.0 mg/dL			
Glu 4.44 mmol/L	tHb 152 g/L	Creat 0.44 mmol/L			
Lac 2.0 mmol/L	tHb 9.4 mmol/L	Creat 442 µmol/L			
Lac 17.8 mg/dL					
BUN 10 mg/dL					
BUN 3.6 mmol/L					
Creat 1.0 mg/dL					
Creat 0.09 mmol/L					
Creat 88.4 µmol/L					

(EN)

Description

The Cartridge contains the following aqueous solutions within individual bags: 5 Calibrators (Calibrator A, B, C, D, and F), Reference Solution, and a flexible waste bag. Each bag includes a filament with a septum. The exposed bag filaments are arranged in a line along the rear of the cartridge box. The septa are pierced during the insertion of the cartridge into the analyser. A change kit consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation.

Included Use:

For the quantitative determination of Na⁺, K⁺, Cl⁻, iCa, iMg, BUN (Urea), Glucose, Lactate, Creatinine, pH, Po₂, PO₂ %, Hct, Hb, O₂ Hb, COHb, MetHb, HHb, and iBil in animal whole blood using this Stat Profile Prime Plus VET Analyzer.

Methodology

Refer to Stat Profile Prime Plus VET Analyzer Instructions For Use Manual for Methodologies and Principles.

Composition

Calibrator bag solution volumes are listed in this table. Contains no constituents of human origin. However, good laboratory practices should be followed during handling of these materials. (REF: NCCLS DOCUMENT M29-T2).

Warnings and Cautions

Refer to Stat Profile Prime Plus VET Analyzer Instructions For Use Manual for complete directions. Follow standard practices required for handling laboratory reagents. Discard Cartridge in accordance with local guidelines.

Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

Storage
Store at 15 - 30°C. Expiration date is printed on the Calibrator Cartridge.

Directions For Use

Activation of Cartridge

Cal A and Cal C must be charged with creatinine syringe prior to installation of the Calibrator Cartridge. The two syringes are labeled and color coded to correspond to the filaments on the carton for Cal A and Cal C.

- Hold the syringe with lip side down and remove protective cap.
- Attach one of the preloaded needle assemblies to the syringe. Remove the protective cap from the needle.
- Match the color and label of the syringe to the appropriate filament and insert needle.
- Slowly depress syringe plunger until the contents are dispensed.
- DO NOT PULL BACK ON THE PLUNGER TO FLUSH CONTENTS OF SYRINGE.
- Remove needle/syringe assembly from filament and discard in an appropriate sharps container.
- Repeat Steps 1-5 for the next Cal.
- Mix Cartridge well by gently shaking for 1 minute.

Refer to the Stat Profile Prime Plus VET Analyzer Instructions For Use Manual for complete directions. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use.

Traceability of Standards

Analytics traced to NIST Standard Reference Materials wherever possible.

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