# Stat Profile Prime Plus® Calibrator Cartridge With Creatinine

Stat Profile Prime Plus® Kalibrierkassette mit Kreatinin, Stat Profile Prime Plus® Φυσιγγίο μέσου βαθμονόμησης με κρεατινίνη, Cartucho de calibrador con creatinina Stat Profile Prime Plus®, Cartouche de calibrateurs Stat Profile Prime Plus® avec créatinine, Cartuccia calibratore con creatinina Stat Profile Prime Plus®, Cartucho de calibrador Stat Profile Prime Plus® com creatinina, Stat Profile Prime Plus® kalibrátor patron kreatininel, סט אפריליזר כבוי עם מומנטני תוצאות Stat Profile Prime Plus®, Stat Profile Prime Plus® キャリブレーターカートリッジ (ケラチニン), 크레이티닌 사용 Stat Profile Prime Plus® 보정기 카트리지, Stat Profile Prime Plus® 含肌酐的校准试剂盒

LOT 23047023

2024-07-26

CAL A	≥ 495 mL	CAL B	≥ 250 mL	CAL C	≥ 245 mL	CAL D	≥ 250 mL	CAL F	≥ 850 mL	SOLN R	≥ 500 mL
HCO <sub>3</sub>	11 mmol/L	pH	6.834 @ 37°C	HCO <sub>3</sub>	22 mmol/L	pH	7.354 @ 37°C	Po <sub>2</sub>	105.5 mmHg	KCl	2 mol/L
Na <sup>+</sup>	135.5 mmol/L	H <sup>+</sup>	146.55 nmol/L	Cl <sup>-</sup>	69.0 mmol/L	H <sup>+</sup>	44.26 nmol/L	Po <sub>2</sub>	14.03 kPa		<span style="border: 1px solid black; padding: 0 2px;">WASTE CONTAINER</span> W
K <sup>+</sup>	4.0 mmol/L	Na <sup>+</sup>	77.0 mmol/L	Glu	200 mg/dL	Lac	10.0 mmol/L				
Cl <sup>-</sup>	108.0 mmol/L	K <sup>+</sup>	10.0 mmol/L	Glu	11.10 mmol/L	Lac	89.0 mg/dL				
iCa	1.04 mmol/L	iCa	1.92 mmol/L	BUN	50 mg/dL						
iCa	4.16 mg/dL	iCa	7.68 mg/dL	BUN	17.9 mmol/L						
iMg	0.50 mmol/L	iMg	1.50 mmol/L	Creat	5.0 mg/dL						
iMg	1.22 mg/dL	iMg	3.63 mg/dL	Creat	0.44 mmol/L						
Glu	80 mg/dL	tHb	15.2 g/dL	Creat	442 µmol/L						
Glu	4.44 mmol/L	tHb	152 g/L								
Lac	2.0 mmol/L	tHb	9.4 mmol/L								
Lac	17.8 mg/dL										
BUN	10 mg/dL										
BUN	3.6 mmol/L										
Creat	1.0 mg/dL										
Creat	0.09 mmol/L										
Creat	88.4 µmol/L										

**EN**

**Description**  
The Cartridge contains the following aqueous solutions within individual bags: 5 Calibrators (Calibrator A, B, C, D, and F), Reference Solution, and a flexible waste bag. Each bag includes a filament in a tube. The exposed bag filaments are arranged in a line along the rear of the cartridge box. The septa are pierced during the insertion of the cartridge into the instrument. A charge kit of three syringes filled with creatinine solution is included for cartridge activation prior to installation.

**Intended Use**  
For the quantitative determination of Na<sup>+</sup>, K<sup>+</sup>, Cl<sup>-</sup>, Ca<sup>2+</sup>, Mg BUN (Urea), Glucose, Lactate, Creatinine, pH, Pco<sub>2</sub>, PO<sub>2</sub>, SO<sub>4</sub><sup>2-</sup>, Hct, Hb, Hbf, O<sub>2</sub>CO<sub>2</sub>, COHb, MetHb, Hrhb, and tBil\* in human blood using the Stat Profile Prime Plus Analyzer.

**Methodology**  
Refer to Stat Profile Prime Plus Analyzer Instructions For Use Manual for Methodologies and Principles.

**Contraindication**  
Calibrator and solution volumes are listed in the table. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF: NCCLS DOCUMENT M29-T2).

**Warnings and Cautions**

Intended for In vitro diagnostic use.

Refer to the Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory reagents. Discard Cartridge in accordance with local guidelines.

**Waste Bag Contains Preservative (ethylene, skin, and respiratory irritant).**

**Storage**

Store at 15 - 30 °C. Expiration date is printed on the Calibrator Cartridge.

**Directions For Use****Activation of Cartridge**

Cal A and Cal C must be charged with creatinine syringe prior to installation of the Calibrator Cartridge. The two syringes are labeled and coded to correspond to the filaments on the cartridge. Cal A and C.

1. Hold the syringe with tip side down and remove protective cap.
2. Attach one of the enclosed needle assemblies to the syringe. Remove the protective cover from the needle.
3. Insert the color-coded end of the syringe to the appropriate filament and insert needle.
4. Slowly draw the syringe plunger until the contents are dispensed.  
DO NOT PUSH BACK ON THE PLUNGER TO FLUSH CONTENTS OF SYRINGE.
5. Remove needle/syringe assembly from filament and discard in an appropriate sharps container.
6. Repeat Steps 1-5 for the next Cal.
7. Mix Cartridge well by gently shaking for 1 minute.

Refer to the Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use.

\*Not available in the USA or for Point-of-Care/Near-Patient Testing use.

**Traceability of Standards**  
Analytes traced to NIST Standard Reference Materials wherever possible.

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Nova Biomedical GmbH Hessenring 13A, Geb. G 64546 Mörfelden-Walldorf Germany
Waltham, MA 02454 U.S.A.