Stat Profile Prime Plus® Calibrator Cartridge With Creatinine

Stat Profile Prime Plus® Kalibrerkassette mit Kreatinin, Stat Profile Prime Plus® Φυσιγγίο μέσου βαθμονόμησης με κρεατινίνη, Cartucho de calibrador con creatinina Stat Profile Prime Plus®, Cartouche de calibrateurs Stat Profile Prime Plus® avec créatinine, Cartuccia calibratore con creatinina Stat Profile Prime Plus®, Cartucho de calibrador Stat Profile Prime Plus® com creatinina, Stat Profile Prime Plus® kalibrátor patron kreatininnel, קריינין סטט פראילס Stat Profile Prime Plus®, Stat Profile Prime Plus® キャリブレータカートリッジ (クリアチニン), 크레아티닌 사용 Stat Profile Prime Plus® 보정기 카트리지, Stat Profile Prime Plus® 含肌酐的校准试剂盒

LOT

23068042

2024-08-13

| CAL A | $\geq 495 \text{ mL}$ | CAL B | $\geq 250 \text{ mL}$ | CAL C | $\geq 245 \text{ mL}$ | CAL D | $\geq 250 \text{ mL}$ | CAL F | $\geq 850 \text{ mL}$ | SOLN R | $\geq 500 \text{ mL}$ |
|------------------|------------------------|---------------|-----------------------|------------------|-----------------------|--------------|-----------------------|---------------|-----------------------|-----------------|-----------------------|
| HCO_3^- | 11 mmol/L | pH | 6.834 @ 37°C | HCO_3^- | 22 mmol/L | pH | 7.354 @ 37°C | PO_2 | 105.7 mmHg | KCl | 2 mol/L |
| Na^+ | 135.5 mmol/L | H^+ | 146.55 nmol/L | Cl^- | 69.0 mmol/L | H^+ | 44.26 nmol/L | PO_2 | 14.06 kPa | WASTE CONTAINER | W |
| K^+ | 4.0 mmol/L | Na^+ | 77.0 mmol/L | Glu | 200 mg/dL | | | | | | |
| Cl^- | 108.0 mmol/L | K^+ | 10.0 mmol/L | Glu | 11.10 mmol/L | | | | | | |
| iCa | 1.04 mmol/L | iCa | 1.92 mmol/L | Lac | 10.0 mmol/L | | | | | | |
| iCa | 4.16 mg/dL | iCa | 7.68 mg/dL | Lac | 89.0 mg/dL | | | | | | |
| iMg | 0.50 mmol/L | iMg | 1.50 mmol/L | BUN | 50 mg/dL | | | | | | |
| iMg | 1.22 mg/dL | iMg | 3.63 mg/dL | BUN | 17.9 mmol/L | | | | | | |
| Glu | 80 mg/dL | tHb | 15.2 g/dL | Creat | 5.0 mg/dL | | | | | | |
| Glu | 4.44 mmol/L | tHb | 152 g/L | Creat | 0.44 mmol/L | | | | | | |
| Lac | 2.0 mmol/L | tHb | 9.4 mmol/L | Creat | 442 $\mu\text{mol/L}$ | | | | | | |
| Lac | 17.8 mg/dL | | | | | | | | | | |
| BUN | 10 mg/dL | | | | | | | | | | |
| BUN | 3.6 mmol/L | | | | | | | | | | |
| Creat | 1.0 mg/dL | | | | | | | | | | |
| Creat | 0.09 mmol/L | | | | | | | | | | |
| Creat | 88.4 $\mu\text{mol/L}$ | | | | | | | | | | |

(EN)

Description

The Cartridge contains the following aqueous solutions within individual bags: 5 Calibrators (Calibrators A, B, C, D, and F), Reference Solution, and a flexible waste bag. Each bag includes a needle assembly. The fittings are arranged in a line along the rear of the cartridge body. The septa are pierced during the insertion of the cartridge into the analyzer. A charge kit consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation.

Intended Use

For the quantitative determination of Na^+ , K^+ , Cl^- , iCa, iMg, BUN (Urea), Glucose, Lactate, Creatinine, pH, PCO_2 , PO_2 , SO_4^{2-} , Hct , tHb, HbF, O₂Hb, COHb, MetHb, HbF, and Bil^+ in human blood using the Stat Profile Prime Plus Analyzer.

Methodology

Refer to Stat Profile Prime Plus Analyzer Instructions for Use Manual for Methodologies and Principles.

Composition

Calibrator and solution volumes are listed in the table. Contains no constituents of human origin, however, good laboratory practice should be followed during handling of these materials. (REF: NCCLS DOCUMENT M29-T2).

Warnings and Cautions

Intended for *in vitro* diagnostic use.

Refer to the Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions. Follow standard practices required for handling laboratory reagents. Discard Cartridge in accordance with local guidelines.

Waste Bag Contains Preservative (eye, skin, and respiratory irritant).

Storage

Store at 15 - 30°C. Expiration date is printed on the Calibrator Cartridge.

Directions For Use

Activation of Cartridge

Cal A and Cal C must be charged with creatinine syringe prior to installation of the Calibrator Cartridge. The two syringes are labeled and color coded to correspond to the fittings on the cartridge. Cal A and Cal C:

- Hold the syringe with tip down and remove protective cap.
- Attach one of the enclosed needle assemblies to the syringe. Remove the protective cover from the needle.
- Match the color and label of the syringes to the appropriate fitting and insert needle.
- DO NOT PULL BACK ON THE PLUNGER TO FLUSH CONTENTS OF SYRINGE.
- Remove needle/syringe assembly from fitting and discard in an appropriate sharps container.
- Repeat Steps 1-5 for the next Cal.
- Mix Cartridge well by gently shaking for 1 minute.

Refer to the Stat Profile Prime Plus Analyzer Instructions for Use Manual for complete directions. Calibrator Cartridge must be inverted gently several times to thoroughly mix contents. Cartridge is ready for use.

*Not available in the USA or for Point-of-Care/Near-Patient Testing use.

Traceability of Standards

Analytes traced to NIST Standard Reference Materials wherever possible.

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Nova Biomedical GmbH Hessenring 13A, Geb. G 64546 Mörfelden-Walldorf Germany
Waltham, MA 02454 U.S.A.