

Attachment 3



**Certificate of Conformance & Analysis**

Nova Biomedical certifies that the following Nova Vet Product has been built and tested in compliance with ISO 13485:2016 and Medical Device Single Audit Program (MDSAP) and meets Nova Biomedical product specifications.

Description: TEST STRIP 10OCT CA NOVAMAX  
 Part Number: 43955  
 BOM Rev.: C  
 Lot Number: 1020222292  
 Expiration Date: 2024-10-19  
 Date of Manufacture: 2022-10-19  
 UDI: (01)20385480 439552  
 (11) 221019  
 (17) 241019  
 (10) 1020222292  
 Quantity: 280  
 Date of Release: 2023-06-05

Release Test Results Summary					
Target Level (mg/dL)	YSI Mean (mg/dL)	Mean (mg/dL)	SD (mg/dL) /CV (%)	SD/CV Specification	P/F
40-80	56.6	61	2.7	≤ 9 mg/dL	P
120-180	174.0	169	3.2	Fail If Both CV of this level and combined CV (L2-L5) >5.0%	P
220-300	297.8	290	3.6	Fail If Both CV of this level and combined CV (L2-L5) >5.0%	P
320-440	426.8	420	4.0	Fail If Both CV of this level and combined CV (L2-L5) >5.0%	P
450-550	488.5	478	3.8	Fail If Both CV of this level and combined CV (L2-L5) >5.0%	P
Total within ISO Zone and %:			(per (ISO15197:2013))	Total within ISO Zone ≥95%	P
<u>360/100%</u>					

**NOTE:** Test data in table above taken from completed and approved Data Sheets from TP 139-501 for this lot number.

Mar Brooke  
 Quality Control Manager Signature  
 Title: In Quality Mgr  
 Date: 6/5/2023

Nova Biomedical Corporation 200 Prospect Street Waltham, MA 02454