

OI-922

Page 62 of 66

Rev. BH

Vialed and Packaged Test Strip Release

Test Strips

nova
biomedical

CERTIFICATE OF COMPLIANCE

This product has been manufactured in compliance with ISO 13485:2016 and Medical Device Single Audit Program (MDSAP) standards and meets Nova Biomedical product specifications.

Customer:	NIA																		
Description:	TEST STRIPS, NOVA VET KETONE, 25 CT., PACKAGED																		
Catalog Number:	52587																		
Nova Biomedical P/N	52587																		
<table><thead><tr><th>Target Range</th><th>MPE</th><th>SD/CV (%)</th></tr></thead><tbody><tr><td>0.1-0.8 mmol/L</td><td>+/- 0.3 mmol/L</td><td>≤ 0.3 mmol/L</td></tr><tr><td>1.5-2.5 mmol/L</td><td>+/- 20%</td><td>≤ 15%</td></tr><tr><td>3.5-5.0 mmol/L</td><td>+/- 20%</td><td>≤ 15%</td></tr></tbody></table>	Target Range	MPE	SD/CV (%)	0.1-0.8 mmol/L	+/- 0.3 mmol/L	≤ 0.3 mmol/L	1.5-2.5 mmol/L	+/- 20%	≤ 15%	3.5-5.0 mmol/L	+/- 20%	≤ 15%	<table><tbody><tr><td>PASS</td><td>FAIL</td></tr><tr><td>PASS</td><td>FAIL</td></tr><tr><td>PASS</td><td>FAIL</td></tr></tbody></table>	PASS	FAIL	PASS	FAIL	PASS	FAIL
Target Range	MPE	SD/CV (%)																	
0.1-0.8 mmol/L	+/- 0.3 mmol/L	≤ 0.3 mmol/L																	
1.5-2.5 mmol/L	+/- 20%	≤ 15%																	
3.5-5.0 mmol/L	+/- 20%	≤ 15%																	
PASS	FAIL																		
PASS	FAIL																		
PASS	FAIL																		
Qty:	518																		
Lot No.	1760223195																		
Mfg. Date:	2023-07-14																		
Expiration Date:	2025-07-14																		
Unique Device Identifier (UDI):	(01)20385480 525873 (11) 230714 (17) 250714 (10) 1760223195																		

Marc Broome
Quality Control Manager

7/28/2023
Date

Nova Biomedical Corporation

200 Prospect Street

Waltham, MA 02454