








5

Assay Data Sheet - Expected Ranges  
BioProfile® FLEX Controls - LEVEL 5

5

Lot Number:  325748		Exp. Date: 2025-03	
Parameter	Units	Expected Range	
		Lower Limit	Upper Limit
Gln	mmol/L	 4.43	 5.93
Glu	mmol/L	 1.78	 2.38
Gluc	g/L	 8.07	 12.07
Lac	g/L	 4.48	 5.98



200 Prospect Street • Waltham, MA U.S.A. 02454-9141 • (781) 864-0800