

Attachment 3



**Certificate of Conformance & Analysis**

Nova Biomedical certifies that the following Nova Vet Product has been built and tested in compliance with ISO 13485:2016 and Medical Device Single Audit Program (MDSAP) and meets Nova Biomedical product specifications.

Description: TEST STRIP 25CT NOVA VET GLU  
 Part Number: 53069  
 BOM Rev.: C  
 Lot Number: 1750424024  
 Expiration Date: 2026-01-24  
 Date of Manufacture: 2024-01-24  
 UDI: (01)20385480 530693  
 (11) 240124  
 (17) 260124  
 (10) 1750424024  
 Quantity: 468  
 Date of Release: 2024-02-08

Release Test Results Summary					
Target Level (mg/dL)	YSI Mean (mg/dL)	Mean (mg/dL)	SD (mg/dL) /CV (%)	SD/CV Specification	P/F
50-100	94.3	111	4.3	≤ 5.0 mg/dL	P
100-140	130.8	138	3.9	≤5.5%	P
160-200	189.8	184	3.0	≤5.0%	P
260-300	283.0	260	3.7	≤5.0%	P
360-400	380.0	342	3.1	≤5.0%	P
Combined MPE (%): <u>2.8%</u>				Specification: <10%	P

**NOTE:** Test data in table above taken from completed and approved Data Sheets from TP 152-058 for this lot number.

Marie Bourne

Quality Control Manager Signature

Title: QC Manager

Date: 3/1/24

Nova Biomedical Corporation 200 Prospect Street Waltham, MA 02454