



**CERTIFICATE OF ANALYSIS**

BioPro Flex Deproteinizing Sol'n Cartridge  
Lot Number 22180004

Ref Number 45192  
Exp Date 2023-06

|         | <u>Spec</u>  | <u>Pass/Fail</u> |
|---------|--------------|------------------|
| NaOCl % | 0.70 - 0.90  | Pass             |
| pH      | 9.00 - 11.00 | Pass             |

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U.S.A in accordance with ISO 13485:2016 Medical Devices Quality Management Systems Requirements, Medical Device Single Audit Program (MDSAP), and conforms to the indicated test specifications. The following analytes are traceable to NIST SRM Materials: pH

Note: Acceptance specifications for this part number are lot dependent and subject to change by the manufacturer

Approval

QC Inspector  
Title

7/29/2022  
Date