



**CERTIFICATE OF ANALYSIS**

Nova StatStrip Glucose and  $\beta$ -Ketone Control Solution Level 1  
Lot Number 0421032301

Ref Number 46947  
Exp Date 2023-08-01  
Date of Manufacture 2021-02-01  
UDI (01)00385480469470(11)210201(17)230801(10)0421032301

	Spec	Mean Result
Glucose mg/dL	57-65	63
Ketone mmol/L	0.1-0.2	0.1

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U.S.A. in accordance with ISO 13485:2016 Medical Devices Quality Management Systems Requirements, Medical Device Single Audit Program (MDSAP), and conforms to the indicated test specifications.

Note: Acceptance specifications for this part number are lot dependent and subject to change by the manufacturer

Approval

QC Inspector  
Title

5/7/2021  
Date