

OI-922

Page 37 of 66

Rev. BL

Vialed and Packaged Test Strip Release

Test Strips

nova
biomedical

CERTIFICATE OF COMPLIANCE

This product has been manufactured in compliance with ISO 13485:2016 and Medical Device Single Audit Program (MDSAP) standards and meets Nova Biomedical product specifications.

Customer:	N/A
Description:	TEST STRIPS, STATSTRIP CAPILLARY GLUCOSE ROW, PACKAGED 100 COUNT
Catalog Number:	61500
Nova Biomedical P/N	61500
Target Range (mg/dL) SD(mg/dL) /CV (%)	
30-50 ≤ 4 mg/dL	PASS / FAIL
51-110 ≤ 6.0%	PASS / FAIL
111-150 ≤ 5.0%	PASS / FAIL
151-250 ≤ 5.0%	PASS / FAIL
251-400 ≤ 5.0%	PASS / FAIL
Accuracy Criteria 1: Total within +/-12 mg/dL or +/- 12% of reference method ≥ 95%	PASS / FAIL
Accuracy Criteria 2: Total within +/-15 mg/dL or +/- 15% of reference method ≥ 98%	PASS / FAIL
Qty:	4,997
Lot No.	0324323309
Mfg. Date:	2024-11-18
Expiration Date:	2027-05-18
Unique Device Identifier (UDI):	(01)20385480615000 (11)241118 (17)270518 (10)0324323309

Mark Broome
Quality Control Manager

12/10/24
Date

Nova Biomedical Corporation

200 Prospect Street Waltham, MA 02454