

Attachment 3



**Certificate of Conformance & Analysis**

Nova Biomedical certifies that the following Nova Vet Product has been built and tested in compliance with ISO 13485:2016 and Medical Device Single Audit Program (MDSAP) and meets Nova Biomedical product specifications.

Description: Test Strip 25ct Novavet Glu  
 Part Number: 53069  
 BOM Rev.: C  
 Lot Number: 1750225009  
 Expiration Date: 2027-01-09  
 Date of Manufacture: 2025-01-09  
 UDI: (01)20385480 530693  
 (11) 250109  
 (17) 270109  
 (10) 1750225009  
 Quantity: 270  
 Date of Release: 2025-01-20

Release Test Results Summary					
Target Level (mg/dL)	YSI Mean (mg/dL)	Mean (mg/dL)	SD (mg/dL) / CV (%)	SD/CV Specification	P/F
50-100	99	97	3.6	≤ 5.0 mg/dL	P
100-140	130	131	2.7	≤ 5.5%	P
160-200	175	184	3.1	≤ 5.0%	P
260-300	249	264	2.9	≤ 5.0%	P
360-400	351	360	2.9	≤ 5.0%	P
Combined MPE (%):	4.1%			Specification: < 10%	P

**NOTE:** Test data in table above taken from completed and approved Data Sheets from TP 152-058 for this lot number.

\_\_\_\_\_  
 Quality Control Manager Signature  
 Title: Sr. Mfg. Manager  
 Date: 1/20/25

Nova Biomedical Corporation 200 Prospect Street Waltham, MA 02454