



CERTIFICATE OF COMPLIANCE

This product has been manufactured in compliance with ISO 13485:2016 and Medical Device Single Audit Program (MDSAP) standards and meets Nova Biomedical product specifications.

Customer:	N/A																
Description:	TEST STRIPS, NOVA VET KETONE, 25 CT., PACKAGED																
Catalog Number:	52587																
Nova Biomedical P/N	52587																
<table border="1"> <thead> <tr> <th>Target Range</th> <th>MPE</th> <th>SD/CV (%)</th> <th></th> </tr> </thead> <tbody> <tr> <td>0.1-0.8 mmol/L</td> <td>+/- 0.3 mmol/L</td> <td>≤ 0.1 mmol/L</td> <td><u>PASS</u> FAIL</td> </tr> <tr> <td>1.5-2.5 mmol/L</td> <td>+/- 20%</td> <td>≤ 10%</td> <td><u>PASS</u> FAIL</td> </tr> <tr> <td>3.5-5.0 mmol/L</td> <td>+/- 20%</td> <td>≤ 10%</td> <td><u>PASS</u> FAIL</td> </tr> </tbody> </table>	Target Range	MPE	SD/CV (%)		0.1-0.8 mmol/L	+/- 0.3 mmol/L	≤ 0.1 mmol/L	<u>PASS</u> FAIL	1.5-2.5 mmol/L	+/- 20%	≤ 10%	<u>PASS</u> FAIL	3.5-5.0 mmol/L	+/- 20%	≤ 10%	<u>PASS</u> FAIL	
Target Range	MPE	SD/CV (%)															
0.1-0.8 mmol/L	+/- 0.3 mmol/L	≤ 0.1 mmol/L	<u>PASS</u> FAIL														
1.5-2.5 mmol/L	+/- 20%	≤ 10%	<u>PASS</u> FAIL														
3.5-5.0 mmol/L	+/- 20%	≤ 10%	<u>PASS</u> FAIL														
Qty:	504																
Lot No.	1760225073																
Mfg. Date:	2025-03-14																
Expiration Date:	2027-03-14																
Unique Device Identifier (UDI):	(01)20385480 <u>525873</u> (11) <u>250314</u> (17) <u>270314</u> (10) <u>1760225073</u>																

Marc Broome 3/31/25
 Quality Control Manager Date

Nova Biomedical Corporation 200 Prospect Street Waltham, MA 02454