



CERTIFICATE OF COMPLIANCE

This product has been manufactured in compliance with ISO 13485:2016 and Medical Device Single Audit Program (MDSAP) standards and meets Nova Biomedical product specifications.

Customer:	NIA	
Description:	TEST STRIPS, STATSTRIP A GLUCOSE INT'L, PACKAGED 100 COUNT	
Catalog Number:	59352	
Nova Biomedical P/N	59352	
Target Range (mg/dL)	SD(mg/dL) /CV (%)	
30-50	≤ 4 mg/dL	PASS / FAIL
51-110	≤ 6.0%	PASS / FAIL
111-150	≤ 5.0%	PASS / FAIL
151-250	≤ 5.0%	PASS / FAIL
251-400	≤ 5.0%	PASS / FAIL
Accuracy Criteria 1: Total within +/-12 mg/dL or +/- 12% of reference method ≥ 95%		PASS / FAIL
Accuracy Criteria 2: Total within +/-15 mg/dL or +/- 15% of reference method ≥ 98%		PASS / FAIL
Qty:	958	
Lot No.	1125086309	
Mfg. Date:	2025-03-27	
Expiration Date:	2027-09-27	
Unique Device Identifier (UDI):	(01)20385480593520 (11)250327 (17)270927 (10)1125086309	

Man Boone
 Quality Control Manager

4/8/25
 Date

Nova Biomedical Corporation

200 Prospect Street

Waltham, MA 02454