



CERTIFICATE OF COMPLIANCE

This product has been manufactured in compliance with ISO 13485:2016 and Medical Device Single Audit Program (MDSAP) standards and meets Nova Biomedical product specifications.

Customer:	N/A	
Description:	TEST STRIPS, STATSTRIP™ GLUCOSE, PACKAGED 1800 COUNT	
Catalog Number:	42214	
Nova Biomedical P/N	55890	
Target Range (mg/dL)	SD(mg/dL) /CV (%)	
30-50	≤ 4 mg/dL	(PASS) / FAIL
51-110	≤ 6.0%	(PASS) / FAIL
111-150	≤ 5.0%	(PASS) / FAIL
151-250	≤ 5.0%	(PASS) / FAIL
251-400	≤ 5.0%	(PASS) / FAIL
Accuracy Criteria 1: Total within +/-12 mg/dL or +/- 12% of reference method ≥ 95%		(PASS) / FAIL
Accuracy Criteria 2: Total within +/-15 mg/dL or +/- 15% of reference method ≥ 98%		(PASS) / FAIL
Qty:	4,689	
Lot No.	0325113249	
Mfg. Date:	2025-04-23	
Expiration Date:	2027-04-23	
Unique Device Identifier (UDI):	(01)20385480422141 (11) 2 5 0 4 2 3 (17) 2 7 0 4 2 3 (10) 0 3 2 5 1 1 3 2 4 9	

Quality Control Manager

5/2/25
Date

Nova Biomedical Corporation

200 Prospect Street Waltham, MA 02454