

Attachment 3



Certificate of Conformance & Analysis

Nova Biomedical certifies that the following Nova Vet Product has been built and tested in compliance with ISO 13485:2016 and Medical Device Single Audit Program (MDSAP) and meets Nova Biomedical product specifications.

Description: 25ct Nova Vet Glu
 Part Number: 53069
 BOM Rev.: C
 Lot Number: 1750225174
 Expiration Date: 2027-06-23
 Date of Manufacture: 2025-06-23
 UDI: (01)20385480 530693
 (11) 250623
 (17) 270623
 (10) 1750225174
 Quantity: 566
 Date of Release: 2025-07-17

Release Test Results Summary					
Target Level (mg/dL)	YSI Mean (mg/dL)	Mean (mg/dL)	SD (mg/dL) /CV (%)	SD/CV Specification	P/F
50-100	88.1	98	2.8	≤ 5.0 mg/dL	P
100-140	126.3	135	4.0	≤5.5%	P
160-200	198.3	199	2.6	≤5.0%	P
260-300	275.5	267	2.6	≤5.0%	P
360-400	366.8	358	2.6	≤5.0%	P
Combined MPE (%):	4.7 %			Specification: <10%	P

NOTE: Test data in table above taken from completed and approved Data Sheets from TP 152-058 for this lot number.

Mark Brown

Quality Control Manager Signature

Title: QC Manager

Date: 7/17/25

Nova Biomedical Corporation 200 Prospect Street Waltham, MA 02454