



CERTIFICATE OF ANALYSIS

BioProfile FLEX 2 CDV Reagent Pack
Lot Number 25317019

Ref Number 55297
Exp Date 2026-10-29
Date of Manufacture 2025-10-29

	CDV Bottle Pack L/N:	25317018	Exp Date:	2026-10-29	
Soln SF L/N	25304018		Soln DP L/N	25303025	
	<u>Spec</u>	<u>Result</u>		<u>Spec</u>	<u>Result</u>
pH	6.900 - 7.100	Pass	Depro Soln	N/A	Present
Cal C L/N	25302052		Soln T L/N	25304017	
	<u>Spec</u>	<u>Result</u>		<u>Spec</u>	<u>Result</u>
3% Detergent	N/A	Present	0.4% Trypan	N/A	Present

CDV Reagent Cartridge L/N: 25317016 Exp Date: 2027-02-05
 Soln PF L/N 25309016
 Spec Pass/Fail
 pH 6.900 - 7.100 Pass

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U.S.A in with ISO 13485:2016 Medical Devices Quality Management Systems Requirements, Medical Device Single Audit Program and conforms to the indicated test specifications. The following analytes are traceable to NIST SRM Materials: pH.

Note: Acceptance specifications for this part number are lot dependent and subject to change by the manufacturer

Approval

QC Inspector
Title

12/5/2025
Date