



**CERTIFICATE OF COMPLIANCE**

This product has been manufactured in compliance with ISO 13485:2016 and Medical Device Single Audit Program (MDSAP) standards and meets Nova Biomedical product specifications.

Customer:	NIPRO DIAGNOSTICS INC.		
Description:	TEST STRIPS, NIPRO STATSTRIP™ KETONE, PACKAGED 20 CT		
Catalog Number:	54288		
Nova Biomedical P/N	54288		
<u>Target Range</u>	<u>MPE</u>	<u>SD/CV (%)</u>	PASS / FAIL PASS / FAIL PASS / FAIL
0.1-0.8 mmol/L	+/- 0.3 mmol/L	≤ 0.1 mmol/L	
1.5-2.5 mmol/L	+/- 20%	≤ 10%	
3.5-5.0 mmol/L	+/- 20%	≤ 10%	
Qty:	834		
Lot No.	2225372249		
Mfg. Date:	2025.11.18		
Expiration Date:	2027.11.18		

*Maureen*  
 Quality Control Manager

12/2/25  
 Date

Nova Biomedical Corporation

200 Prospect Street Waltham, MA 02454