

Attachment 3



**Certificate of Conformance & Analysis**

Nova Biomedical certifies that the following Nova Vet Product has been built and tested in compliance with ISO 13485:2016 and Medical Device Single Audit Program (MDSAP) and meets Nova Biomedical product specifications.

Description: TEST STRIP 25CT NOVA VET GIU  
 Part Number: 53069  
 BOM Rev.: C  
 Lot Number: 1750226051  
 Expiration Date: 2028-02-20  
 Date of Manufacture: 2026-02-20  
 UDI: (01)20385480 530693  
 (11) 260220  
 (17) 280220  
 (10) 1750226051  
 Quantity: 210  
 Date of Release: 2026-03-11

Release Test Results Summary

Target Level (mg/dL)	YSI Mean (mg/dL)	Mean (mg/dL)	SD (mg/dL) /CV (%)	SD/CV Specification	P/F
50-100	107	103	3.1	≤ 5.0 mg/dL	P
100-140	141	139	3.2	≤5.5%	P
160-200	179	181	3.0	≤5.0%	P
260-300	262	262	3.7	≤5.0%	P
360-400	357	377	3.0	≤5.0%	P
Combined MPE (%): <u>3.4</u>				Specification: <10%	P

**NOTE:** Test data in table above taken from completed and approved Data Sheets from TP 152-058 for this lot number.

[Signature]  
 Quality Control Manager Signature  
 Title: QC Manager  
 Date: 3/11/26

Nova Biomedical Corporation 200 Prospect Street Waltham, MA 02454